



# 一孔双钉系列

(圆孔 & 结合孔)

## 锁定接骨板固定系统

(Round-holes & Combi-holes)

LOCKING PLATE FIXATION SYSTEM



江苏双羊医疗器械有限公司

Jiangsu Shuangyang Medical Instrument Co., Ltd.



# Company Introduction

## 企业简介



江苏双羊医疗器械有限公司，创建于2001年，占地面积18000平方米，建筑面积15000多平方米，注册资金2000万。是一家从事骨科植入产品的研发、制造、销售和服务的民族企业，获得多项国家级专利。

公司选择钛及钛合金为原材料，严把质量关，选择国内外知名企业宝钛、ZAPP等作为原材料供应商，并配备世界一流的加工中心、纵切车床、数控铣床、超声波清洗等生产设备和装置，以及万能试验机、电子扭转机、数字投影仪等精密的量具。公司拥有完善的管理体系，先后通过了ISO9001:2015质量管理体系认证，ISO13485:2016医疗器械质量管理体系认证，TUV的CE认证，并于2007年率先通过国家局组织的《医疗器械生产质量管理规范植入性医疗器械实施细则（试点）》检查。

公司针对人体各个骨骼部位的特点，在著名骨科专家、教授、临床医师的悉心指导与帮助下，推出众多主导产品：锁定接骨板固定系统、钛质接骨板固定系统、钛质空心接骨螺钉及垫片、钛质胸肋骨系统、锁定颌面内固定系统、颌面内固定系统、钛质捆绑系统、解剖型钛网系统、基础工具系列等，并拥有专业的配套手术器械包，以满足临床不同的需求。公司的产品因设计可靠、加工精细、医生临床使用方便、患者愈合期短而受到了许多临床医师和患者的高度评价。

中国梦，双羊梦！双羊作为一个有使命感、有责任心、有抱负、有人文情怀的公司的初衷不会变，始终以“以人为本，诚信至上，不断创新，追求卓越”的理念不会变，争创医疗器械领域民族品牌领头羊的决心不会变。愿更多有志者与双羊共同创造美好的未来！



Jiangsu Shuangyang Medical Instrument Co., Ltd. was founded in 2001, covers an area of 18,000 m<sup>2</sup>, including a floor area of over 15,000 m<sup>2</sup>. Its registered capital reaches 20 million Yuan. As a national enterprise dedicated to the R&D, manufacturing, sales and service of orthopaedic implants, we have obtained several national patents.

Titanium and titanium alloys are our raw materials. We perform strict quality control, and select domestic and international famous brands, such as Baotai and ZAPP, as our raw material suppliers. Meanwhile, we are equipped with world class production equipment and devices including machining center, slitting lathe, CNC milling machine, and ultrasonic cleaner, etc., as well as precise measuring implements including

universal tester, electronic torsion tester and digital projector, etc. Thanks to a sophisticated management system, we have acquired ISO9001: 2015 Certificate of Quality Management System, ISO13485:2016 Certificate of Quality Management System for Medical Devices, and CE certificate of TUV. We are also the first to pass the inspection according to the *Enforcement Regulation (Pilot) for Implantable Medical Devices of Good Manufacturing Practice for Medical Devices* organized by the National Bureau in 2007.

Thanks to meticulous guidance and supports from distinguished orthopaedic specialists, professors and clinicians, we have launched numerous leading products customized for different human skeletal parts, including locking bone plate fixation system, titanium bone plate fixation system, titanium cannulated bone screw & gasket, titanium sternocostal system, locking maxillofacial internal fixation system, maxillofacial internal fixation system, titanium binding system, anatomic titanium mesh system, and basic tool series, etc. We also have professional supporting surgical instrument sets to meet various clinical needs. Extensive acclaims have been received from clinicians and patients for our easy-to-use products with reliable design and fine machining, which can bring a short healing period.

China dream and Shuangyang dream! We will stick to our original intention to be a mission-driven, responsible, ambitious and humanistic company, and adhere to our idea of "people orientation, integrity, innovation, and excellence". We are determined to be a leading national brand in the medical Instrument industry. At Shuangyang, we always welcome aspiring talents to co-create a bright future with us.





# Advanced Equipment

## 精良装备

专业的技术支持、先进的设备、严格的工艺规程、  
配套零部件的精心选用、一丝不苟的敬业精神、充分保  
证了产品的质量。

Professional technical support, advanced  
equipment, strict technological procedures,  
carefully selected auxiliary parts, and meticulous  
conscientiousness, all these factors help to ensure  
the quality of our products.





纵切车削中心 The Turning Center

# Quality Control System

## 品控体系



我们坚守稳定的产品质量，恪守产品的表里如一，信奉制造的一丝不苟，从设计、制造、检测，到管理，每一道工序，每一个环节，专家式的管理，皆表露出 ISO9001:2000 的规范化和标准化。

We try our best to ensure reliable product quality and meticulous manufacturing. From design, manufacturing, detection to management, we carry out professional control in each step and each process according to ISO9001:2000 norms and standards.



①



②

1. 原材料管控  
Raw Material Control
2. 设备刀具辅料管控  
Equipment, Cutter & Accessory Control



③



④



⑤

3. 过程能力管控  
Process Capability Control
4. 工装模具管控  
Tooling Control
5. 质量能力管控  
Quality Capacity Control

# Contents

## 目录

### 12-14 锁定 2.0 系列接骨板 2.0 series locking plate

- 14 锁定桡骨小头接骨板  
Radial head locking plate



### 15-17 锁定 2.4 系列接骨板 2.4 series locking plate

- 17 锁定桡骨小头接骨板  
Radial head locking plate



### 18-28 锁定 3.0 系列接骨板 3.0 series locking plate

- 18 锁定鹰嘴接骨板  
Olecranon locking plate
- 19 锁定掌背侧接骨板  
Volar dorsal locking plate
- 20 锁定掌侧接骨板  
Volar locking plate
- 23 锁定 3.0 系列接骨板  
3.0 series locking plate
- 25 锁定 1/3 管型接骨板  
One-third tubular locking plate
- 26 锁定腓骨远端接骨板  
Distal fibular locking plate
- 28 锁定跟骨接骨板  
Calcaneal locking plate





## 30-55 锁定 4.0 系列接骨板

### 4.0 series locking plate

- 30 锁定 4.0 直形接骨板  
4.0 straight locking plate
- 32 锁定骨盆重建接骨板  
Pelvic reconstruction locking plate
- 34 锁定锁骨重建型接骨板（中远端）  
Clavicle reconstruction locking plate (Middle & Distal)
- 35 锁定锁骨重建型接骨板（远端）  
Clavicle reconstruction locking plate (Distal)
- 37 锁定锁骨钩接骨板  
Clavicle hook locking plate
- 38 锁定肱骨颈接骨板  
Neck of humeral locking plate
- 39 万向锁定肱骨颈接骨板  
Multi-axial neck of humeral locking plate
- 40 锁定肱骨下段 Y 形接骨板（前侧）  
Anterior humeral Y shape locking plate
- 41 锁定肱骨下段 Y 形接骨板（后侧）  
Posterior humeral Y shape locking plate
- 42 锁定肱骨远端亚髁接骨板  
Distal humeral sub-condyle locking plate
- 43 锁定肱骨远端内侧接骨板  
Distal medial humeral locking plate
- 44 锁定肱骨远端外侧接骨板  
Distal lateral humeral locking plate
- 46 锁定胫骨平台外侧接骨板  
Lateral tibia plateau locking plate
- 47 万向锁定胫骨平台外侧接骨板  
Multi-axial lateral tibia plateau locking plate
- 48 锁定胫骨平台内侧接骨板  
Medial tibia plateau locking plate
- 49 万向锁定胫骨平台内侧接骨板  
Multi-axial medial tibia plateau locking plate
- 50 锁定胫骨平台正后侧接骨板  
Posterior tibia plateau locking plate

- 52 锁定胫骨平台后内侧接骨板  
Posteromedial tibia plateau locking plate
- 53 锁定胫骨远端后外侧接骨板  
Distal posterior lateral tibia locking plate
- 54 锁定胫骨远端外侧 L 型接骨板  
Distal lateral tibia L shape locking plate
- 55 锁定胫骨远端内侧接骨板  
Distal medial tibia locking plate

## 56-69 锁定 5.0 系列接骨板

### 5.0 series locking plate

- 56 锁定 5.0 窄型接骨板  
5.0 straight narrow locking plate
- 57 锁定 5.0 宽型接骨板  
5.0 straight broad locking plate
- 58 万向锁定股骨颈接骨板 I 型  
Multi-axial neck of femur locking plate ( I type)
- 59 万向锁定股骨颈接骨板 II 型  
Multi-axial neck of femur locking plate ( II type)
- 60 锁定股骨近端接骨板  
Proximal femur locking plate (Combi-holes)
- 61 万向锁定股骨近端接骨板  
Multi-axial proximal femur locking plate
- 62 锁定股骨远端接骨板  
Distal femur locking plate
- 63 万向锁定股骨远端接骨板  
Multi-axial distal femur locking plate
- 64 锁定胫骨平台内侧接骨板  
Medial tibia plateau locking plate



- 65 万向锁定胫骨平台内侧接骨板  
Multi-axial medial tibia plateau locking plate
- 66 锁定胫骨平台外侧接骨板  
Lateral tibia plateau locking plate
- 67 万向锁定胫骨平台外侧接骨板  
Multi-axial lateral tibia plateau locking plate
- 68 锁定胫骨远端外侧接骨板  
Distal lateral tibia locking plate
- 69 万向锁定胫骨远端外侧接骨板  
Multi-axial distal lateral tibia locking plate

- 70-74 钛质锁定接骨螺钉      Titanium locking screw
- 75-79 钛质皮质骨接骨螺钉      Titanium cortex screw
- 80-85 钛质松质骨接骨螺钉      Titanium cancellous screw
- 86-98 钛质空心接骨螺钉      Titanium cannulated screw



# 锁定 2.0 系列接骨板

## 2.0 series locking plate

**适用部位：**指骨、掌骨、趾骨、跖骨、桡骨头

Suitable for phalange, metacarpal, phalanx, metatarsal, head of radius

**配用：**Φ2.0 锁定钉，Φ2.0 皮质骨螺钉

**Matching** Φ2.0 locking screw, Φ2.0 cortex screw

### 直形 Straight

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.06.02.04011010   | 直形               | 4 孔 | 24mm |
| ★ 10.06.02.05011010 | 直形               | 5 孔 | 30mm |
| 10.06.02.06011010   | 直形               | 6 孔 | 36mm |
| 10.06.02.07011010   | 直形               | 7 孔 | 42mm |



### 90° L 形 90° L-plate

| 订货编码 Order code.     | 规格 Specification |     |      |
|----------------------|------------------|-----|------|
| 10.06.02.041115010   | 90° L 形 左 (L)    | 4 孔 | 20mm |
| 10.06.02.042115010   | 90° L 形 右 (R)    | 4 孔 | 20mm |
| ★ 10.06.02.051115010 | 90° L 形 左 (L)    | 5 孔 | 26mm |
| 10.06.02.052115010   | 90° L 形 右 (R)    | 5 孔 | 26mm |



## 110° L 形 110° L-plate

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.06.02.04116010   | 110° L 形 左 (L)   | 4 孔 | 20mm |
| 10.06.02.04216010   | 110° L 形 右 (R)   | 4 孔 | 20mm |
| ★ 10.06.02.05116010 | 110° L 形 左 (L)   | 5 孔 | 26mm |
| 10.06.02.05216010   | 110° L 形 右 (R)   | 5 孔 | 26mm |



## T 形 T-plate

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.06.02.04019010   | T 形              | 4 孔 | 20mm |
| 10.06.02.05019010   | T 形              | 5 孔 | 26mm |
| 10.06.02.06019010   | T 形              | 6 孔 | 32mm |
| ★ 10.06.02.07019010 | T 形              | 7 孔 | 38mm |



## T 形 T-plate

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 00.06.02.05019015   | T 形              | 5 孔 | 23mm |
| 00.06.02.06019015   | T 形              | 6 孔 | 29mm |
| 00.06.02.07019015   | T 形              | 7 孔 | 35mm |
| ★ 00.06.02.08019015 | T 形              | 8 孔 | 41mm |





## 锁定桡骨小头接骨板

Radial head locking plate

**适用部位：**桡骨小头

Suitable for capitulum radius

**配用：**Φ2.0 锁定钉，Φ2.0 皮质骨螺钉

**Matching** Φ2.0 locking screw, Φ2.0 cortex screw



| 订货编码 Order code.    | 规格 Specification |      |      |
|---------------------|------------------|------|------|
| 10.06.02.07019000   | T 形              | 7 孔  | 29mm |
| * 10.06.02.08019000 | T 形              | 8 孔  | 35mm |
| 10.06.02.09019000   | T 形              | 9 孔  | 41mm |
| 10.06.02.10019000   | T 形              | 10 孔 | 47mm |

# 锁定 2.4 系列接骨板

## 2.4 series locking plate

**适用部位：**指骨、掌骨、跖骨、小儿尺骨及桡骨

Suitable for phalange, metacarpal, metatarsal, pediatric ulna and radius

**配用：**Φ2.4 锁定钉，Φ2.4 皮质骨螺钉

**Macthing** Φ2.4 locking screw, Φ2.4 cortex screw

### 直形 Straight

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.14.03.04011015   | 直形               | 4 孔 | 36mm |
| 10.14.03.05011015   | 直形               | 5 孔 | 45mm |
| ★ 10.14.03.06011015 | 直形               | 6 孔 | 54mm |
| 10.14.03.08011015   | 直形               | 8 孔 | 72mm |





## 90° L 形 90° L-plate

| 订货编码 Order code.    | 规格 Specification |     |      |  |
|---------------------|------------------|-----|------|--|
| 10.14.03.04115015   | 90° L 型左 (L)     | 4 孔 | 27mm |  |
| 10.14.03.04215015   | 90° L 型右 (R)     | 4 孔 | 27mm |  |
| ★ 10.14.03.05115015 | 90° L 型左 (L)     | 5 孔 | 36mm |  |
| 10.14.03.05215015   | 90° L 型右 (R)     | 5 孔 | 36mm |  |



## 110° L 形 110° L-plate

| 订货编码 Order code.    | 规格 Specification |     |      |  |
|---------------------|------------------|-----|------|--|
| 10.14.03.04116015   | 110° L 型左 (L)    | 4 孔 | 30mm |  |
| 10.14.03.04216015   | 110° L 型右 (R)    | 4 孔 | 30mm |  |
| ★ 10.14.03.05116015 | 110° L 型左 (L)    | 5 孔 | 39mm |  |
| 10.14.03.05216015   | 110° L 型右 (R)    | 5 孔 | 39mm |  |



## T 形 T-plate

| 订货编码 Order code.    | 规格 Specification |      |      |  |
|---------------------|------------------|------|------|--|
| 10.14.03.05019015   | T 形              | 5 孔  | 27mm |  |
| 10.14.03.06019015   | T 形              | 6 孔  | 36mm |  |
| ★ 10.14.03.07019015 | T 形              | 7 孔  | 45mm |  |
| 10.14.03.08019015   | T 形              | 8 孔  | 54mm |  |
| 10.14.03.10019015   | T 形              | 10 孔 | 72mm |  |



# 锁定桡骨小头接骨板

## Radial head locking plate



**适用部位：**桡骨小头

Suitable for capitulum radius

**配用：**Φ2.4 锁定钉，Φ2.4 皮质骨螺钉

**Matching** Φ2.4 locking screw, Φ2.4 cortex screw

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.14.03.05019020   | T 形              | 5 孔 | 37mm |
| ★ 10.14.03.06019020 | T 形              | 6 孔 | 46mm |
| 10.14.03.07019020   | T 形              | 7 孔 | 55mm |



## 锁定 3.0 系列接骨板

### 3.0 series locking plate

#### 锁定鹰嘴接骨板

#### Olecranon locking plate

**适用部位：**尺骨鹰嘴

Suitable for ulna olecranon

**配用：**Φ3.0 皮质骨螺钉，Φ3.0 锁定钉

**Matching** Φ3.0 cortex screw, Φ3.0 locking screw



| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.16.04102000   | 左 L              | 4 孔 | 88mm  |
| 10.14.16.04202000   | 右 R              | 4 孔 | 88mm  |
| * 10.14.16.05102000 | 左 L              | 5 孔 | 103mm |
| 10.14.16.05202000   | 右 R              | 5 孔 | 103mm |
| 10.14.16.06102000   | 左 L              | 6 孔 | 115mm |
| 10.14.16.06202000   | 右 R              | 6 孔 | 115mm |
| 10.14.16.07102000   | 左 L              | 7 孔 | 127mm |
| 10.14.16.07202000   | 右 R              | 7 孔 | 127mm |
| 10.14.16.08102000   | 左 L              | 8 孔 | 139mm |
| 10.14.16.08202000   | 右 R              | 8 孔 | 139mm |

## 锁定掌背侧接骨板

Volar dorsal locking plate



**适用部位:** 桡骨远端掌背侧 Suitable for distal volar dorsal radius

**配用:**  $\Phi 3.0$  皮质骨螺钉,  $\Phi 3.0$  锁定钉

**Matching**  $\Phi 3.0$  cortex screw,  $\Phi 3.0$  locking screw

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.14.18.03102000   | 左 L              | 3 孔 | 51mm |
| 10.14.18.03202000   | 右 R              | 3 孔 | 51mm |
| 10.14.18.04102000   | 左 L              | 4 孔 | 63mm |
| 10.14.18.04202000   | 右 R              | 4 孔 | 63mm |
| * 10.14.18.05102000 | 左 L              | 5 孔 | 75mm |
| 10.14.18.05202000   | 右 R              | 5 孔 | 75mm |
| 10.14.18.06102000   | 左 L              | 6 孔 | 87mm |
| 10.14.18.06202000   | 右 R              | 6 孔 | 87mm |



## 锁定掌侧接骨板 Volar locking plate

**适用部位:** 桡骨远端掌侧  
Suitable for distal volar radius

### 斜头型 Oblique head type

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.20.03104000   | 左 L              | 3 孔 | 57mm  |
| 10.14.20.03204000   | 右 R              | 3 孔 | 57mm  |
| 10.14.20.04104000   | 左 L              | 4 孔 | 69mm  |
| 10.14.20.04204000   | 右 R              | 4 孔 | 69mm  |
| ★ 10.14.20.05104000 | 左 L              | 5 孔 | 81mm  |
| 10.14.20.05204000   | 右 R              | 5 孔 | 81mm  |
| 10.14.20.06104000   | 左 L              | 6 孔 | 93mm  |
| 10.14.20.06204000   | 右 R              | 6 孔 | 93mm  |
| 10.14.20.08104000   | 左 L              | 8 孔 | 117mm |
| 10.14.20.08204000   | 右 R              | 8 孔 | 117mm |



### 远端型 Distal type

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.11.21.03102077   | 左 L              | 3 孔 | 47mm |
| 10.11.21.03202077   | 右 R              | 3 孔 | 47mm |
| 10.11.21.04102077   | 左 L              | 4 孔 | 58mm |
| 10.11.21.04202077   | 右 R              | 4 孔 | 58mm |
| ★ 10.11.21.05102077 | 左 L              | 5 孔 | 69mm |
| 10.11.21.05202077   | 右 R              | 5 孔 | 69mm |
| 10.11.21.06102077   | 左 L              | 6 孔 | 80mm |
| 10.11.21.06202077   | 右 R              | 6 孔 | 80mm |



配用:  $\Phi$ 3.0 皮质骨螺钉,  $\Phi$ 3.0 锁定钉

Matching  $\Phi$ 3.0 cortical screw,  $\Phi$ 3.0 locking screw

### 斜头型 ( 小号 )

### Oblique head type (small)

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| ★ 10.14.20.03105011 | 左 L              | 3 孔 | 64mm  |
| 10.14.20.03205011   | 右 R              | 3 孔 | 64mm  |
| 10.14.20.04105011   | 左 L              | 4 孔 | 79mm  |
| 10.14.20.04205011   | 右 R              | 4 孔 | 79mm  |
| 10.14.20.05105011   | 左 L              | 5 孔 | 91mm  |
| 10.14.20.05205011   | 右 R              | 5 孔 | 91mm  |
| 10.14.20.06105011   | 左 L              | 6 孔 | 103mm |
| 10.14.20.06205011   | 右 R              | 6 孔 | 103mm |



### 斜头型 ( 大号 )

### Oblique head type (large)

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| ★ 10.14.20.03105021 | 左 L              | 3 孔 | 64mm  |
| 10.14.20.03205021   | 右 R              | 3 孔 | 64mm  |
| 10.14.20.04105021   | 左 L              | 4 孔 | 79mm  |
| 10.14.20.04205021   | 右 R              | 4 孔 | 79mm  |
| 10.14.20.05105021   | 左 L              | 5 孔 | 91mm  |
| 10.14.20.05205021   | 右 R              | 5 孔 | 91mm  |
| 10.14.20.06105021   | 左 L              | 6 孔 | 103mm |
| 10.14.20.06205021   | 右 R              | 6 孔 | 103mm |





梅花型 (小号) Torx type (small)

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.20.03105112   | 左 L              | 3 孔 | 64mm  |
| 10.14.20.03205112   | 右 R              | 3 孔 | 64mm  |
| ★ 10.14.20.04105112 | 左 L              | 4 孔 | 79mm  |
| 10.14.20.04205112   | 右 R              | 4 孔 | 79mm  |
| 10.14.20.06105112   | 左 L              | 6 孔 | 103mm |
| 10.14.20.06205112   | 右 R              | 6 孔 | 103mm |



梅花型 (大号) Torx type (large)

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.20.03105122   | 左 L              | 3 孔 | 64mm  |
| 10.14.20.03205122   | 右 R              | 3 孔 | 64mm  |
| ★ 10.14.20.04105122 | 左 L              | 4 孔 | 79mm  |
| 10.14.20.04205122   | 右 R              | 4 孔 | 79mm  |
| 10.14.20.06105122   | 左 L              | 6 孔 | 103mm |
| 10.14.20.06205122   | 右 R              | 6 孔 | 103mm |



## 锁定 3.0 系列接骨板

### 3.0 series locking plate

**适用部位:** 尺骨、桡骨

Suitable for ulna and radius

**配用:**  $\Phi$ 3.0 皮质骨螺钉,  $\Phi$ 3.0 锁定钉

Matching  $\Phi$ 3.0 cortical screw,  $\Phi$ 3.0 locking screw

#### 90° L 形 90° L-plate

| 订货编码 Order code.    | 规格 Specification |       |      |  |
|---------------------|------------------|-------|------|--|
| 10.14.06.03115001   | 90° L 形 (L)      | 2-3 孔 | 51mm |  |
| 10.14.06.03215001   | 90° L 形 (R)      | 2-3 孔 | 51mm |  |
| 10.14.06.04115001   | 90° L 形 (L)      | 2-4 孔 | 63mm |  |
| 10.14.06.04215001   | 90° L 形 (R)      | 2-4 孔 | 63mm |  |
| 10.14.06.03115002   | 90° L 形 (L)      | 3-3 孔 | 51mm |  |
| 10.14.06.03215002   | 90° L 形 (R)      | 3-3 孔 | 51mm |  |
| ★ 10.14.06.04115002 | 90° L 形 (L)      | 3-4 孔 | 63mm |  |
| 10.14.06.04215002   | 90° L 形 (R)      | 3-4 孔 | 63mm |  |



#### 110° L 形 110° L-plate

| 订货编码 Order code.    | 规格 Specification |       |      |  |
|---------------------|------------------|-------|------|--|
| 10.14.06.03116001   | 110° L 形 (L)     | 2-3 孔 | 54mm |  |
| 10.14.06.03216001   | 110° L 形 (R)     | 2-3 孔 | 54mm |  |
| 10.14.06.04116001   | 110° L 形 (L)     | 2-4 孔 | 66mm |  |
| 10.14.06.04216001   | 110° L 形 (R)     | 2-4 孔 | 66mm |  |
| 10.14.06.03116002   | 110° L 形 (L)     | 3-3 孔 | 57mm |  |
| 10.14.06.03216002   | 110° L 形 (R)     | 3-3 孔 | 57mm |  |
| ★ 10.14.06.04116002 | 110° L 形 (L)     | 3-4 孔 | 69mm |  |
| 10.14.06.04216002   | 110° L 形 (R)     | 3-4 孔 | 69mm |  |





### T 形 T-plate

| 订货编码 Order code.    | 规格 Specification |       |      |
|---------------------|------------------|-------|------|
| 10.14.06.03019001   | T 形              | 2-3 孔 | 51mm |
| 10.14.06.04019001   | T 形              | 2-4 孔 | 63mm |
| 10.14.06.05019001   | T 形              | 2-5 孔 | 75mm |
| 10.14.06.06019001   | T 形              | 2-6 孔 | 87mm |
| 10.14.06.03019002   | T 形              | 3-3 孔 | 52mm |
| ★ 10.14.06.04019002 | T 形              | 3-4 孔 | 64mm |
| 10.14.06.05019002   | T 形              | 3-5 孔 | 76mm |
| 10.14.06.06019002   | T 形              | 3-6 孔 | 88mm |



### 桡骨远端外侧 Distal lateral radius

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| 10.11.06.05011000   | 直形               | 5 孔 | 48mm |
| ★ 10.11.06.06011000 | 直形               | 6 孔 | 57mm |

### 尺桡骨远端 Distal ulna and radius

| 订货编码 Order code.    | 规格 Specification |      |       | 宽度 width |
|---------------------|------------------|------|-------|----------|
| 10.14.06.04011025   | 直形               | 4 孔  | 47mm  | 8mm      |
| ★ 10.14.06.06011025 | 直形               | 6 孔  | 71mm  |          |
| 10.14.06.08011025   | 直形               | 8 孔  | 95mm  |          |
| 10.14.06.10011025   | 直形               | 10 孔 | 119mm |          |
| 10.14.06.12011025   | 直形               | 12 孔 | 143mm |          |



## 尺桡骨远端 Distal ulna and radius

| 订货编码 Order code.    | 规格 Specification |      |       | 宽度 width |
|---------------------|------------------|------|-------|----------|
| 10.14.06.06011125   | 直形               | 6 孔  | 77mm  | 9mm      |
| ★ 10.14.06.07011125 | 直形               | 7 孔  | 89mm  |          |
| 10.14.06.08011125   | 直形               | 8 孔  | 101mm |          |
| 10.14.06.10011125   | 直形               | 10 孔 | 125mm |          |
| 10.14.06.12011125   | 直形               | 12 孔 | 149mm |          |



## 锁定 1/3 管型接骨板

## One-third tubular locking plate

适用部位：腓骨

Suitable for fibula

配用：

Φ3.0 皮质骨螺钉

Φ3.0 锁定钉

Matching:

Φ3.0 cortex screw

Φ3.0 locking screw

| 订货编码 Order code.    | 规格 Specification |       |
|---------------------|------------------|-------|
| 10.11.17.05000000   | 5 孔              | 61mm  |
| ★ 10.11.17.06000000 | 6 孔              | 73mm  |
| 10.11.17.07000000   | 7 孔              | 85mm  |
| 10.11.17.08000000   | 8 孔              | 97mm  |
| 10.11.17.09000000   | 9 孔              | 109mm |
| 10.11.17.10000000   | 10 孔             | 121mm |





## 锁定腓骨远端接骨板

Distal fibular locking plate

**适用部位：**腓骨骨端

Suitable for capitulum fibula

**配用：**

Φ3.0 皮质骨螺钉

Φ3.0 锁定钉

**Matching:**

Φ3.0 cortex screw

Φ3.0 locking screw



**I 型 - 前外侧 Anterior lateral**

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.35.04101000   | 左 L              | 4 孔 | 85mm  |
| 10.14.35.04201000   | 右 R              | 4 孔 | 85mm  |
| ★ 10.14.35.05101000 | 左 L              | 5 孔 | 98mm  |
| 10.14.35.05201000   | 右 R              | 5 孔 | 98mm  |
| 10.14.35.06101000   | 左 L              | 6 孔 | 111mm |
| 10.14.35.06201000   | 右 R              | 6 孔 | 111mm |
| 10.14.35.07101000   | 左 L              | 7 孔 | 124mm |
| 10.14.35.07201000   | 右 R              | 7 孔 | 124mm |
| 10.14.35.08101000   | 左 L              | 8 孔 | 137mm |
| 10.14.35.08201000   | 右 R              | 8 孔 | 137mm |

## II型 - 后外侧 Posterior lateral

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.35.04102000   | 左 L              | 4 孔 | 83mm  |
| 10.14.35.04202000   | 右 R              | 4 孔 | 83mm  |
| * 10.14.35.05102000 | 左 L              | 5 孔 | 95mm  |
| 10.14.35.05202000   | 右 R              | 5 孔 | 95mm  |
| 10.14.35.06102000   | 左 L              | 6 孔 | 107mm |
| 10.14.35.06202000   | 右 R              | 6 孔 | 107mm |
| 10.14.35.08102000   | 左 L              | 8 孔 | 131mm |
| 10.14.35.08202000   | 右 R              | 8 孔 | 131mm |



## III型 - 外侧 Lateral

| 订货编码 Order code.    | 规格 Specification |       |
|---------------------|------------------|-------|
| 10.14.35.04003000   | 4 孔              | 79mm  |
| * 10.14.35.05003000 | 5 孔              | 91mm  |
| 10.14.35.06003000   | 6 孔              | 103mm |
| 10.14.35.08003000   | 8 孔              | 127mm |





## 锁定跟骨接骨板 Calcaneal locking plate

**适用部位：**跟骨

Suitable for calcaneus

**配用：**

- Φ4.0 松质骨螺钉
- Φ3.0 锁定钉
- Φ3.0 皮质骨螺钉



**Matching:**

- Φ4.0 cancellous screw
- Φ3.0 locking screw
- Φ3.0 cortex screw

**I 型 - 爪型 claw type**

| 订货编码 Order code.    | 规格 Specification |      |      |
|---------------------|------------------|------|------|
| ★ 10.11.36.12129260 | 左 L              | 12 孔 | 60mm |
| 10.11.36.12229260   | 右 R              | 12 孔 | 60mm |
| 10.11.36.14129269   | 左 L              | 14 孔 | 69mm |
| 10.11.36.14229269   | 右 R              | 14 孔 | 69mm |

## 锁定跟骨接骨板

### Calcaneal locking plate



#### II型 - 环抱型 loop type

| 订货编码 Order code.    | 规格 Specification |      |      |
|---------------------|------------------|------|------|
| 10.11.36.10129200   | 左 L              | 10 孔 | 48mm |
| 10.11.36.10229200   | 右 R              | 10 孔 | 48mm |
| ★ 10.11.36.12129200 | 左 L              | 12 孔 | 58mm |
| 10.11.36.12229200   | 右 R              | 12 孔 | 58mm |
| 10.11.36.14129200   | 左 L              | 14 孔 | 65mm |
| 10.11.36.14229200   | 右 R              | 14 孔 | 65mm |



#### 平板型 Flat type

| 订货编码 Order code.    | 规格 Specification |      |      |
|---------------------|------------------|------|------|
| 10.11.36.12102260   | 左                | 12 孔 | 60mm |
| 10.11.36.12202260   | 右                | 12 孔 | 60mm |
| ★ 10.11.36.12102269 | 左                | 12 孔 | 69mm |
| 10.11.36.12202269   | 右                | 12 孔 | 69mm |



## 锁定 4.0 系列接骨板

4.0 series locking plate

### 锁定 4.0 直形接骨板

4.0 straight locking plate

**适用部位：**尺骨、桡骨、肱骨、小儿胫骨或股骨

Suitable for ulna, radius, humerus, pediatric tibia or femur

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.5 cortex screw

Φ4.0 locking screw



**宽型 Broad type**

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.07.06011000   | 直形               | 6 孔  | 87mm  |
| 10.14.07.07011000   | 直形               | 7 孔  | 100mm |
| ★ 10.14.07.08011000 | 直形               | 8 孔  | 113mm |
| 10.14.07.09011000   | 直形               | 9 孔  | 126mm |
| 10.14.07.10011000   | 直形               | 10 孔 | 139mm |
| 10.14.07.12011000   | 直形               | 12 孔 | 165mm |



**适用部位：**尺骨、桡骨、肱骨、小儿胫骨或股骨

Suitable for ulna, radius, humerus, pediatric tibia or femur

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.5 cortex screw

Φ4.0 locking screw



**窄型 Narrow type**

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.07.06011001   | 直形               | 6 孔  | 87mm  |
| 10.14.07.07011001   | 直形               | 7 孔  | 100mm |
| * 10.14.07.08011001 | 直形               | 8 孔  | 113mm |
| 10.14.07.09011001   | 直形               | 9 孔  | 126mm |
| 10.14.07.10011001   | 直形               | 10 孔 | 139mm |
| 10.14.07.12011001   | 直形               | 12 孔 | 165mm |



## 锁定骨盆重建接骨板

Pelvic reconstruction locking plate

**适用部位：**髌臼、骨盆

Suitable for acetabulum, pelvic

### 重建直形 Straight reconstruction



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.11.07.06030000   | 重建直形             | 6 孔  | 72mm  |
| 10.11.07.07030000   | 重建直形             | 7 孔  | 84mm  |
| ★ 10.11.07.08030000 | 重建直形             | 8 孔  | 96mm  |
| 10.11.07.09030000   | 重建直形             | 9 孔  | 108mm |
| 10.11.07.10030000   | 重建直形             | 10 孔 | 120mm |
| 10.11.07.12030000   | 重建直形             | 12 孔 | 144mm |



### 重建弧形 Arc reconstruction

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| ★ 10.11.07.06031000 | 重建弧形             | 6 孔  | 69mm  |
| 10.11.07.08031000   | 重建弧形             | 8 孔  | 89mm  |
| 10.11.07.10031000   | 重建弧形             | 10 孔 | 107mm |
| 10.11.07.12031000   | 重建弧形             | 12 孔 | 122mm |
| 10.11.07.14031000   | 重建弧形             | 14 孔 | 133mm |
| 10.11.07.16031000   | 重建弧形             | 16 孔 | 139mm |

配用：

Φ4.0 松质骨螺钉

Φ4.0 锁定钉

Matching:

Φ4.0 cancellous screw

Φ4.0 locking screw

## 重建直形 II Straight reconstruction II



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.11.07.06030002   | 重建直形             | 6 孔  | 72mm  |
| 10.11.07.07030002   | 重建直形             | 7 孔  | 84mm  |
| * 10.11.07.08030002 | 重建直形             | 8 孔  | 96mm  |
| 10.11.07.09030002   | 重建直形             | 9 孔  | 108mm |
| 10.11.07.10030002   | 重建直形             | 10 孔 | 120mm |
| 10.11.07.12030002   | 重建直形             | 12 孔 | 144mm |



## 重建弧形 II Arc reconstruction II

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.11.07.06031002   | 重建弧形             | 6 孔  | 71mm  |
| * 10.11.07.08031002 | 重建弧形             | 8 孔  | 92mm  |
| 10.11.07.10031002   | 重建弧形             | 10 孔 | 110mm |
| 10.11.07.12031002   | 重建弧形             | 12 孔 | 124mm |
| 10.11.07.14031002   | 重建弧形             | 14 孔 | 134mm |
| 10.11.07.16031002   | 重建弧形             | 16 孔 | 140mm |



## 锁定锁骨重建型接骨板（中远端）

Clavicle reconstruction locking plate (middle & Distal)

**适用部位：**锁骨中远端

Suitable for middle and distal clavicle

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 松质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.5 cortex screw

Φ4.0 cancellous screw

Φ4.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.10.06128022   | 中段左 L            | 6 孔  | 86mm  |
| 10.14.10.06228022   | 中段右 R            | 6 孔  | 86mm  |
| 10.14.10.07128022   | 中段左 L            | 7 孔  | 100mm |
| 10.14.10.07228022   | 中段右 R            | 7 孔  | 100mm |
| 10.14.10.08128022   | 中段左 L            | 8 孔  | 112mm |
| 10.14.10.08228022   | 中段右 R            | 8 孔  | 112mm |
| 10.14.10.10128022   | 中段左 L            | 10 孔 | 135mm |
| 10.14.10.10228022   | 中段右 R            | 10 孔 | 135mm |
| ★ 10.14.10.06128021 | 远端左 L            | 6 孔  | 83mm  |
| 10.14.10.06228021   | 远端右 R            | 6 孔  | 83mm  |
| 10.14.10.07128021   | 远端左 L            | 7 孔  | 96mm  |
| 10.14.10.07228021   | 远端右 R            | 7 孔  | 96mm  |
| 10.14.10.08128021   | 远端左 L            | 8 孔  | 109mm |
| 10.14.10.08228021   | 远端右 R            | 8 孔  | 109mm |
| 10.14.10.10128021   | 远端左 L            | 10 孔 | 134mm |
| 10.14.10.10228021   | 远端右 R            | 10 孔 | 134mm |

## 锁定锁骨重建型接骨板（远端）

Clavicle reconstruction locking plate (Distal)



**适用部位：**锁骨远端

Suitable for distal clavicle

**配用：**

Φ3.0 皮质骨螺钉  
 Φ3.0 锁定钉  
 Φ3.5 皮质骨螺钉  
 Φ4.0 松质骨螺钉  
 Φ4.0 锁定钉

**Matching:**

Φ3.0 cortex screw  
 Φ3.0 locking screw  
 Φ3.5 cortex screw  
 Φ4.0 cancellous screw  
 Φ4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.10.04128001   | 远端左 L            | 4 孔 | 75mm  |
| 10.14.10.04228001   | 远端右 R            | 4 孔 | 75mm  |
| ★ 10.14.10.06128001 | 远端左 L            | 6 孔 | 102mm |
| 10.14.10.06228001   | 远端右 R            | 6 孔 | 102mm |
| 10.14.10.07128001   | 远端左 L            | 7 孔 | 115mm |
| 10.14.10.07228001   | 远端右 R            | 7 孔 | 115mm |



## 锁定锁骨重建型接骨板（远端）

Clavicle reconstruction locking plate (Distal)



**适用部位：**锁骨远端

Suitable for distal clavicle

**配用：**

- Φ3.0 皮质骨螺钉
- Φ3.0 锁定钉
- Φ3.5 皮质骨螺钉
- Φ4.0 松质骨螺钉
- Φ4.0 锁定钉

**Matching:**

- Φ3.0 cortex screw
- Φ3.0 locking screw
- Φ3.5 cortex screw
- Φ4.0 cancellous screw
- Φ4.0 locking screw

**冠状解剖型 Corona anatomical**

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.11.10.04128001   | 远端左 L            | 4 孔 | 73mm  |
| 10.11.10.04228001   | 远端右 R            | 4 孔 | 73mm  |
| 10.11.10.06128001   | 远端左 L            | 6 孔 | 95mm  |
| ★ 10.11.10.06228001 | 远端右 R            | 6 孔 | 95mm  |
| 10.11.10.07128001   | 远端左 L            | 7 孔 | 106mm |
| 10.11.10.07228001   | 远端右 R            | 7 孔 | 106mm |
| 10.11.10.08128001   | 远端左 L            | 8 孔 | 115mm |
| 10.11.10.08228001   | 远端右 R            | 8 孔 | 115mm |

## 锁定锁骨钩接骨板

### Clavicle hook locking plate



钩高 15.5mm Hook height

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.11.05100100   | 左 L              | 5 孔 | 76mm  |
| 10.14.11.05200100   | 右 R              | 5 孔 | 76mm  |
| ★ 10.14.11.06100100 | 左 L              | 6 孔 | 89mm  |
| 10.14.11.06200100   | 右 R              | 6 孔 | 89mm  |
| 10.14.11.07100100   | 左 L              | 7 孔 | 102mm |
| 10.14.11.07200100   | 右 R              | 7 孔 | 102mm |
| 10.14.11.08100100   | 左 L              | 8 孔 | 115mm |
| 10.14.11.08200100   | 右 R              | 8 孔 | 115mm |

**适用部位：**锁骨肩峰端

Suitable for acromial end clavicle

**配用：**

- Φ3.5 皮质骨螺钉
- Φ4.0 松质骨螺钉
- Φ4.0 锁定钉

**Matching:**

- Φ3.5 cortex screw
- Φ4.0 cancellous screw
- Φ4.0 locking screw

钩高 18mm Hook height

| 订货编码 Order code.  | 规格 Specification |     |       |
|-------------------|------------------|-----|-------|
| 10.14.11.05100118 | 左 L              | 5 孔 | 76mm  |
| 10.14.11.05200118 | 右 R              | 5 孔 | 76mm  |
| 10.14.11.06100118 | 左 L              | 6 孔 | 89mm  |
| 10.14.11.06200118 | 右 R              | 6 孔 | 89mm  |
| 10.14.11.07100118 | 左 L              | 7 孔 | 102mm |
| 10.14.11.07200118 | 右 R              | 7 孔 | 102mm |
| 10.14.11.08100118 | 左 L              | 8 孔 | 115mm |
| 10.14.11.08200118 | 右 R              | 8 孔 | 115mm |



## 锁定肱骨颈接骨板

Neck of humerus locking plate

**适用部位:** 肱骨近端

Suitable for proximal humerus

**配用:**  $\Phi 3.5$  皮质骨螺钉,  $\Phi 4.0$  松质骨螺钉,  $\Phi 4.0$  锁定钉

**Matching:**  $\Phi 3.5$  cortex screw,  $\Phi 4.0$  cancellous screw,  $\Phi 4.0$  locking screw



| 订货编码 Order code.    | 规格 Specification |       |
|---------------------|------------------|-------|
| 10.14.12.03001300   | 3 孔              | 89mm  |
| 10.14.12.04001300   | 4 孔              | 102mm |
| * 10.14.12.05001300 | 5 孔              | 115mm |
| 10.14.12.06001300   | 6 孔              | 128mm |
| 10.14.12.07001300   | 7 孔              | 141mm |
| 10.14.12.08001300   | 8 孔              | 154mm |
| 10.14.12.10001300   | 10 孔             | 180mm |
| 10.14.12.12001300   | 12 孔             | 206mm |

## 万向锁定肱骨颈接骨板

Multi-axial neck of humerus locking plate

**适用部位：**肱骨近端

Suitable for proximal humerus

**配用：**Φ3.5 皮质骨螺钉，Φ4.0 松质骨螺钉，Φ4.0 锁定钉

**Matching：**Φ3.5 cortex screw，Φ4.0 cancellous screw，Φ4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |       |
|---------------------|------------------|-------|
| 10.14.13.03001000   | 3 孔              | 89mm  |
| 10.14.13.04001000   | 4 孔              | 102mm |
| * 10.14.13.05001000 | 5 孔              | 115mm |
| 10.14.13.06001000   | 6 孔              | 128mm |
| 10.14.13.07001000   | 7 孔              | 141mm |
| 10.14.13.08001000   | 8 孔              | 154mm |
| 10.14.13.10001000   | 10 孔             | 180mm |
| 10.14.13.12001000   | 12 孔             | 206mm |





## 锁定肱骨下段 Y 形接骨板 (前侧) Anterior humeral Y-shape locking plate

**适用部位:** 肱骨下段前侧

Suitable for middle lower anterior humerus

**配用:**

- Φ3.5 皮质骨螺钉
- Φ4.0 松质骨螺钉
- Φ4.0 锁定钉

**Matching:**

- Φ3.5 cortex screw
- Φ4.0 cancellous screw
- Φ4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |       |
|---------------------|------------------|-------|
| 10.11.07.06020300   | 6 孔              | 136mm |
| ★ 10.11.07.07020300 | 7 孔              | 152mm |
| 10.11.07.08020300   | 8 孔              | 168mm |
| 10.11.07.09020300   | 9 孔              | 184mm |
| 10.11.07.10020300   | 10 孔             | 200mm |
| 10.11.07.12020300   | 12 孔             | 232mm |



## 锁定肱骨下段 Y 形接骨板 (后侧)

### Posterior humeral Y-shape locking plate

**适用部位:** 肱骨中下段后侧

Suitable for middle lower posterior humerus

**配用:**

- Φ3.5 皮质骨螺钉
- Φ4.0 松质骨螺钉
- Φ4.0 锁定钉

**Matching:**

- Φ3.5 cortex screw
- Φ4.0 cancellous screw
- Φ4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| ★ 10.11.07.07120300 | 左 L              | 7 孔  | 150mm |
| 10.11.07.07220300   | 右 R              | 7 孔  | 150mm |
| 10.11.07.09120300   | 左 L              | 9 孔  | 182mm |
| 10.11.07.09220300   | 右 R              | 9 孔  | 182mm |
| 10.11.07.11120300   | 左 L              | 11 孔 | 214mm |
| 10.11.07.11220300   | 右 R              | 11 孔 | 214mm |





## 锁定肱骨远端亚髁接骨板

Distal humeral sub-condyle locking plate

**适用部位：** 肱骨中下段外侧

Suitable for distal humeral fracture

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 松质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.5 cortex screw

Φ4.0 cancellous screw

Φ4.0 locking screw



| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.14.07.04111000   | 左 L              | 4 孔 | 104mm |
| 10.14.07.04211000   | 右 R              | 4 孔 | 104mm |
| 10.14.07.06111000   | 左 L              | 6 孔 | 139mm |
| 10.14.07.06211000   | 右 R              | 6 孔 | 139mm |
| * 10.14.07.07111000 | 左 L              | 7 孔 | 155mm |
| 10.14.07.07211000   | 右 R              | 7 孔 | 155mm |
| 10.14.07.08111000   | 左 L              | 8 孔 | 174mm |
| 10.14.07.08211000   | 右 R              | 8 孔 | 174mm |
| 10.14.07.09111000   | 左 L              | 9 孔 | 190mm |
| 10.14.07.09211000   | 右 R              | 9 孔 | 190mm |

## 锁定肱骨远端内侧接骨板

Distal medial humerus locking plate

**适用部位：** 肱骨远端内侧

Suitable for distal medial humerus

**配用：**

Φ3.0 皮质骨螺钉  
 Φ3.0 锁定钉  
 Φ3.5 皮质骨螺钉  
 Φ4.0 松质骨螺钉  
 Φ4.0 锁定钉

**Matching:**

Φ3.0 cortex screw  
 Φ3.0 locking screw  
 Φ3.5 cortex screw  
 Φ4.0 cancellous screw  
 Φ4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.15.09100000   | 左 L              | 9 孔  | 103mm |
| 10.14.15.09200000   | 右 R              | 9 孔  | 103mm |
| ★ 10.14.15.11100000 | 左 L              | 11 孔 | 129mm |
| 10.14.15.11200000   | 右 R              | 11 孔 | 129mm |





## 锁定肱骨远端外侧接骨板

Distal lateral humerus locking plate

**适用部位：** 肱骨远端外侧

Suitable for distal lateral humerus

**配用：**

- Φ3.0 皮质骨螺钉
- Φ3.0 锁定钉
- Φ3.5 皮质骨螺钉
- Φ4.0 松质骨螺钉
- Φ4.0 锁定钉


**Matching:**

- Φ3.0 cortex screw
- Φ3.0 locking screw
- Φ3.5 cortex screw
- Φ4.0 cancellous screw
- Φ4.0 locking screw

**I 型 - 外侧支撑型 Lateral support type**

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.14.07100000   | 左 L              | 7 孔  | 123mm |
| 10.14.14.07200000   | 右 R              | 7 孔  | 123mm |
| ★ 10.14.14.09100000 | 左 L              | 9 孔  | 149mm |
| 10.14.14.09200000   | 右 R              | 9 孔  | 149mm |
| 10.14.14.11100000   | 左 L              | 11 孔 | 175mm |
| 10.14.14.11200000   | 右 R              | 11 孔 | 175mm |

| 订货编码 Order code. | 规格 Specification |      |       |
|------------------|------------------|------|-------|
| 99.10.19074.00   | 左 L              | 7 孔  | 72mm  |
| 99.10.19077.00   | 右 R              | 7 孔  | 72mm  |
| 99.10.19075.00   | 左 L              | 9 孔  | 98mm  |
| 99.10.19078.00   | 右 R              | 9 孔  | 98mm  |
| 99.10.19076.00   | 左 L              | 11 孔 | 124mm |
| 99.10.19079.00   | 右 R              | 11 孔 | 124mm |



## 锁定肱骨远端外侧接骨板

Distal lateral humerus locking plate

**适用部位：** 肱骨远端外侧

Suitable for distal lateral humerus

**配用：**

Φ3.0 皮质骨螺钉

Φ3.0 锁定钉

Φ3.5 皮质骨螺钉

Φ4.0 松质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.0 cortex screw

Φ3.0 locking screw

Φ3.5 cortex screw

Φ4.0 cancellous screw

Φ4.0 locking screw

### II 型

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.14.07102000   | 左 L              | 7 孔  | 123mm |
| 10.14.14.07202000   | 右 R              | 7 孔  | 123mm |
| * 10.14.14.09102000 | 左 L              | 9 孔  | 149mm |
| 10.14.14.09202000   | 右 R              | 9 孔  | 149mm |
| 10.14.14.11102000   | 左 L              | 11 孔 | 175mm |
| 10.14.14.11202000   | 右 R              | 11 孔 | 175mm |

| 订货编码 Order code. | 规格 Specification |      |       |
|------------------|------------------|------|-------|
| 99.10.19080.00   | 左 L              | 7 孔  | 72mm  |
| 99.10.19083.00   | 右 R              | 7 孔  | 72mm  |
| 99.10.19081.00   | 左 L              | 9 孔  | 98mm  |
| 99.10.19084.00   | 右 R              | 9 孔  | 98mm  |
| 99.10.19082.00   | 左 L              | 11 孔 | 124mm |
| 99.10.19085.00   | 右 R              | 11 孔 | 124mm |





## 锁定胫骨平台外侧接骨板

Lateral tibia plateau locking plate

**适用部位:** 胫骨平台外侧 Suitable for lateral tibia plateau

**配用:**  $\Phi 3.5$  皮质骨螺钉,  $\Phi 4.0$  松质骨螺钉,  $\Phi 4.0$  锁定钉  
**Matching**  $\Phi 3.5$  cortex screw,  $\Phi 4.0$  cancellous screw,  
 $\Phi 4.0$  locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.30.03103000   | 左 L              | 3 孔  | 86mm  |
| 10.14.30.03203000   | 右 R              | 3 孔  | 86mm  |
| 10.14.30.05103000   | 左 L              | 5 孔  | 118mm |
| 10.14.30.05203000   | 右 R              | 5 孔  | 118mm |
| ★ 10.14.30.07103000 | 左 L              | 7 孔  | 150mm |
| 10.14.30.07203000   | 右 R              | 7 孔  | 150mm |
| 10.14.30.09103000   | 左 L              | 9 孔  | 182mm |
| 10.14.30.09203000   | 右 R              | 9 孔  | 182mm |
| 10.14.30.11103000   | 左 L              | 11 孔 | 214mm |
| 10.14.30.11203000   | 右 R              | 11 孔 | 214mm |
| 10.14.30.13103000   | 左 L              | 13 孔 | 246mm |
| 10.14.30.13203000   | 右 R              | 13 孔 | 246mm |
| 10.14.30.15103000   | 左 L              | 15 孔 | 278mm |
| 10.14.30.15203000   | 右 R              | 15 孔 | 278mm |

## 万向锁定胫骨平台外侧接骨板

### Multi-axial lateral tibia plateau locking plate

**适用部位:** 胫骨平台外侧 Suitable for lateral tibia plateau

**配用:**  $\Phi 3.5$  皮质骨螺钉,  $\Phi 4.0$  松质骨螺钉,  $\Phi 4.0$  锁定钉

**Matching**  $\Phi 3.5$  cortex screw,  $\Phi 4.0$  cancellous screw,  $\Phi 4.0$  locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.31.03103000   | 左 L              | 3 孔  | 86mm  |
| 10.14.31.03203000   | 右 R              | 3 孔  | 86mm  |
| 10.14.31.05103000   | 左 L              | 5 孔  | 118mm |
| 10.14.31.05203000   | 右 R              | 5 孔  | 118mm |
| ★ 10.14.31.07103000 | 左 L              | 7 孔  | 150mm |
| 10.14.31.07203000   | 右 R              | 7 孔  | 150mm |
| 10.14.31.09103000   | 左 L              | 9 孔  | 182mm |
| 10.14.31.09203000   | 右 R              | 9 孔  | 182mm |
| 10.14.31.11103000   | 左 L              | 11 孔 | 214mm |
| 10.14.31.11203000   | 右 R              | 11 孔 | 214mm |
| 10.14.31.13103000   | 左 L              | 13 孔 | 246mm |
| 10.14.31.13203000   | 右 R              | 13 孔 | 246mm |
| 10.14.31.15103000   | 左 L              | 15 孔 | 278mm |
| 10.14.31.15203000   | 右 R              | 15 孔 | 278mm |



## 锁定胫骨平台内侧接骨板 Medial tibia plateau locking plate

**适用部位:** 胫骨平台内侧 Suitable for medial tibia plateau

**配用:**  $\Phi 3.5$  皮质骨螺钉,  $\Phi 4.0$  松质骨螺钉,  $\Phi 4.0$  锁定钉

**Matching:**  $\Phi 3.5$  cortex screw,  
 $\Phi 4.0$  cancellous screw,  $\Phi 4.0$  locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.28.03103000   | 左 L              | 3 孔  | 83mm  |
| 10.14.28.03203000   | 右 R              | 3 孔  | 83mm  |
| 10.14.28.05103000   | 左 L              | 5 孔  | 115mm |
| 10.14.28.05203000   | 右 R              | 5 孔  | 115mm |
| ★ 10.14.28.07103000 | 左 L              | 7 孔  | 147mm |
| 10.14.28.07203000   | 右 R              | 7 孔  | 147mm |
| 10.14.28.09103000   | 左 L              | 9 孔  | 179mm |
| 10.14.28.09203000   | 右 R              | 9 孔  | 179mm |
| 10.14.28.11103000   | 左 L              | 11 孔 | 211mm |
| 10.14.28.11203000   | 右 R              | 11 孔 | 211mm |
| 10.14.28.13103000   | 左 L              | 13 孔 | 243mm |
| 10.14.28.13203000   | 右 R              | 13 孔 | 243mm |
| 10.14.28.15103000   | 左 L              | 15 孔 | 275mm |
| 10.14.28.15203000   | 右 R              | 15 孔 | 275mm |

## 万向锁定胫骨平台内侧接骨板

Multi-axial medial tibia plateau locking plate

**适用部位：** 胫骨平台内侧 Suitable for medial tibia plateau

**配用：**  $\Phi 3.5$  皮质骨螺钉,  $\Phi 4.0$  松质骨螺钉,  $\Phi 4.0$  锁定钉

**Matching:**  $\Phi 3.5$  cortex screw,  
 $\Phi 4.0$  cancellous screw,  $\Phi 4.0$  locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.29.03103000   | 左 L              | 3 孔  | 83mm  |
| 10.14.29.03203000   | 右 R              | 3 孔  | 83mm  |
| 10.14.29.05103000   | 左 L              | 5 孔  | 115mm |
| 10.14.29.05203000   | 右 R              | 5 孔  | 115mm |
| ★ 10.14.29.07103000 | 左 L              | 7 孔  | 147mm |
| 10.14.29.07203000   | 右 R              | 7 孔  | 147mm |
| 10.14.29.09103000   | 左 L              | 9 孔  | 179mm |
| 10.14.29.09203000   | 右 R              | 9 孔  | 179mm |
| 10.14.29.11103000   | 左 L              | 11 孔 | 211mm |
| 10.14.29.11203000   | 右 R              | 11 孔 | 211mm |
| 10.14.29.13103000   | 左 L              | 13 孔 | 243mm |
| 10.14.29.13203000   | 右 R              | 13 孔 | 243mm |
| 10.14.29.15103000   | 左 L              | 15 孔 | 275mm |
| 10.14.29.15203000   | 右 R              | 15 孔 | 275mm |



## 锁定胫骨平台正后侧接骨板 Posterior tibia plateau locking plate

**适用部位：** 胫骨平台正后侧

Suitable for: Posterior tibia plateau

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 松质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.5 cortex screw

Φ4.0 cancellous screw

Φ4.0 locking screw

**I 型**



| 订货编码 Order code.    | 规格 Specification |     |       |  |
|---------------------|------------------|-----|-------|--|
| 10.11.07.03101000   | 左 L              | 3 孔 | 77mm  |  |
| 10.11.07.03201000   | 右 R              | 3 孔 | 77mm  |  |
| 10.11.07.04101000   | 左 L              | 4 孔 | 89mm  |  |
| 10.11.07.04201000   | 右 R              | 4 孔 | 89mm  |  |
| * 10.11.07.05101000 | 左 L              | 5 孔 | 101mm |  |
| 10.11.07.05201000   | 右 R              | 5 孔 | 101mm |  |
| 10.11.07.06101000   | 左 L              | 6 孔 | 113mm |  |
| 10.11.07.06201000   | 右 R              | 6 孔 | 113mm |  |
| 10.11.07.07101000   | 左 L              | 7 孔 | 125mm |  |
| 10.11.07.07201000   | 右 R              | 7 孔 | 125mm |  |

## 锁定胫骨平台正后侧接骨板

### Posterior tibia plateau locking plate

**适用部位：**胫骨平台正后侧

Suitable for: Posterior tibia plateau

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 松质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.5 cortex screw

Φ4.0 cancellous screw

Φ4.0 locking screw



#### II型

| 订货编码 Order code.    | 规格 Specification |     |       |
|---------------------|------------------|-----|-------|
| 10.11.07.03102000   | 左 L              | 3 孔 | 67mm  |
| 10.11.07.03202000   | 右 R              | 3 孔 | 67mm  |
| 10.11.07.04102000   | 左 L              | 4 孔 | 79mm  |
| 10.11.07.04202000   | 右 R              | 4 孔 | 79mm  |
| ★ 10.11.07.05102000 | 左 L              | 5 孔 | 91mm  |
| 10.11.07.05202000   | 右 R              | 5 孔 | 91mm  |
| 10.11.07.06102000   | 左 L              | 6 孔 | 103mm |
| 10.11.07.06202000   | 右 R              | 6 孔 | 103mm |
| 10.11.07.07102000   | 左 L              | 7 孔 | 115mm |
| 10.11.07.07202000   | 右 R              | 7 孔 | 115mm |



## 锁定胫骨平台后内侧接骨板

Posteromedial tibia plateau locking plate

**适用部位：**胫骨平台后内侧

Suitable for: Posteromedial tibia plateau

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 松质骨螺钉

Φ4.0 锁定钉

**Matching:**

Φ3.5 cortex screw

Φ4.0 cancellous screw

Φ4.0 locking screw



| 订货编码 Order code.    | 规格 Specification |       |
|---------------------|------------------|-------|
| 10.14.28.03000000   | 3 孔              | 105mm |
| * 10.14.28.05000000 | 5 孔              | 131mm |
| 10.14.28.07000000   | 7 孔              | 157mm |
| 10.14.28.09000000   | 9 孔              | 185mm |

## 锁定胫骨远端后外侧接骨板

Distal posterolateral tibia locking plate

**适用部位：** 胫骨远端后外侧

Suitable for posterolateral distal tibia

**配用：**

Φ3.5 皮质骨螺钉

Φ4.0 松质骨螺钉

Φ4.0 锁定钉



**Matching:**

Φ3.5 cortex screw

Φ4.0 cancellous screw

Φ4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |     |      |
|---------------------|------------------|-----|------|
| ★ 10.11.07.03120000 | 左 L              | 3 孔 | 62mm |
| 10.11.07.03220000   | 右 R              | 3 孔 | 62mm |
| 10.11.07.05120000   | 左 L              | 5 孔 | 86mm |
| 10.11.07.05220000   | 右 R              | 5 孔 | 86mm |



## 锁定胫骨远端外侧 L 型接骨板

Distal lateral tibia L-shaped locking plate

**适用部位：** 胫骨远端外侧

Suitable for distal lateral tibia

**配用：**

- Φ 3.5 皮质骨螺钉
- Φ 4.0 松质骨螺钉
- Φ 4.0 锁定钉

**Matching**

- Φ 3.5 cortex screw
- Φ 4.0 cancellous screw
- Φ 4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.33.05103000   | 左 L              | 5 孔  | 103mm |
| 10.14.33.05203000   | 右 R              | 5 孔  | 103mm |
| * 10.14.33.07103000 | 左 L              | 7 孔  | 135mm |
| 10.14.33.07203000   | 右 R              | 7 孔  | 135mm |
| 10.14.33.09103000   | 左 L              | 9 孔  | 167mm |
| 10.14.33.09203000   | 右 R              | 9 孔  | 167mm |
| 10.14.33.11103000   | 左 L              | 11 孔 | 199mm |
| 10.14.33.11203000   | 右 R              | 11 孔 | 199mm |
| 10.14.33.13103000   | 左 L              | 13 孔 | 231mm |
| 10.14.33.13203000   | 右 R              | 13 孔 | 231mm |
| 10.14.33.15103000   | 左 L              | 15 孔 | 263mm |
| 10.14.33.15203000   | 右 R              | 15 孔 | 263mm |



## 锁定胫骨远端内侧接骨板

### Distal medial tibia locking plate

**适用部位：** 胫骨远端内侧

Suitable for distal medial tibia

**配用：**

- Φ3.5 皮质骨螺钉
- Φ4.0 松质骨螺钉
- Φ4.0 锁定钉

**Matching**

- Φ3.5 cortex screw
- Φ4.0 cancellous screw
- Φ4.0 locking screw

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.32.05103000   | 左 L              | 5 孔  | 116mm |
| 10.14.32.05203000   | 右 R              | 5 孔  | 116mm |
| ★ 10.14.32.07103000 | 左 L              | 7 孔  | 148mm |
| 10.14.32.07203000   | 右 R              | 7 孔  | 148mm |
| 10.14.32.09103000   | 左 L              | 9 孔  | 180mm |
| 10.14.32.09203000   | 右 R              | 9 孔  | 180mm |
| 10.14.32.11103000   | 左 L              | 11 孔 | 212mm |
| 10.14.32.11203000   | 右 R              | 11 孔 | 212mm |
| 10.14.32.13103000   | 左 L              | 13 孔 | 244mm |
| 10.14.32.13203000   | 右 R              | 13 孔 | 244mm |
| 10.14.32.15103000   | 左 L              | 15 孔 | 276mm |
| 10.14.32.15203000   | 右 R              | 15 孔 | 276mm |





## 锁定 5.0 系列接骨板

5.0 series locking plate

### 锁定 5.0 窄型接骨板

5.0 straight narrow locking plate

**适用部位:** 肱骨、胫骨、小儿股骨

Suitable for humerus, tibia, pediatric femur

**配用:**

Φ4.5 皮质骨螺钉

Φ5.0 锁定钉

**Matching**

Φ4.5 cortex screw

Φ5.0 locking screw

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.08.08011300   | 直形               | 8 孔  | 145mm |
| 10.14.08.09011300   | 直形               | 9 孔  | 161mm |
| 10.14.08.10011300   | 直形               | 10 孔 | 177mm |
| * 10.14.08.12011300 | 直形               | 12 孔 | 209mm |
| 10.14.08.14011300   | 直形               | 14 孔 | 241mm |



## 锁定 5.0 宽型接骨板

### 5.0 straight broad locking plate

**适用部位：**胫骨、股骨

Suitable for tibia, femur

**配用：**

Φ4.5 皮质骨螺钉

Φ5.0 锁定钉

**Matching**

Φ4.5 cortex screw

Φ5.0 locking screw

**解剖型 Anatomical**

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.09.09011301   | 直形               | 9 孔  | 195mm |
| 10.14.09.10011301   | 直形               | 10 孔 | 213mm |
| ★ 10.14.09.12011301 | 直形               | 12 孔 | 249mm |
| 10.14.09.14011301   | 直形               | 14 孔 | 284mm |
| 10.14.09.16011301   | 直形               | 16 孔 | 320mm |





## 万向锁定股骨颈接骨板 I 型

Multi-axial neck of femur locking plate I type

**适用部位：**股骨粗隆骨折

Suitable for femoral trochanter fracture

**配用：**

**Matching**

Φ6.5 空心接骨螺钉

Φ6.5 cannulated screw

Φ6.5 空心锁定钉

Φ6.5 cannulated locking screw

Φ6.5 松质骨螺钉

Φ6.5 cancellous screw

Φ4.5 皮质骨螺钉

Φ4.5 cortex screw

Φ5.0 锁定钉

Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.11.22.03100301   | 左 L              | 3 孔  | 135mm |
| 10.11.22.03200301   | 右 R              | 3 孔  | 135mm |
| 10.11.22.05100301   | 左 L              | 5 孔  | 169mm |
| 10.11.22.05200301   | 右 R              | 5 孔  | 169mm |
| 10.11.22.07100301   | 左 L              | 7 孔  | 203mm |
| 10.11.22.07200301   | 右 R              | 7 孔  | 203mm |
| 10.11.22.09100301   | 左 L              | 9 孔  | 237mm |
| ★ 10.11.22.09200301 | 右 R              | 9 孔  | 237mm |
| 10.11.22.11100301   | 左 L              | 11 孔 | 271mm |
| 10.11.22.11200301   | 右 R              | 11 孔 | 271mm |
| 10.11.22.13100301   | 左 L              | 13 孔 | 304mm |
| 10.11.22.13200301   | 右 R              | 13 孔 | 304mm |

## 万向锁定股骨颈接骨板 II 型

Multi-axial neck of femur locking plate II type

**适用部位：**股骨颈骨折

Suitable for fracture of neck of femur

### 配用：

Φ6.5 空心接骨螺钉

Φ6.5 空心锁定钉

Φ6.5 松质骨螺钉

Φ5.0 锁定钉

### Matching

Φ6.5 cannulated screw

Φ6.5 cannulated locking screw

Φ6.5 cancellous screw

Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification      |     |     |      |
|---------------------|-----------------------|-----|-----|------|
| * 10.11.22.02003000 | 万向<br>Multi-<br>axial | 左右式 | 2 孔 | 66mm |
| 10.11.22.02102000   |                       | 左 L | 2 孔 | 66mm |
| 10.11.22.02202000   |                       | 右 R | 2 孔 | 66mm |



## 锁定股骨近端接骨板

### Proximal femur locking plate

**适用部位：**粗隆下端

Suitable for proximal femur

**配用：**

Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉

**Matching** Φ4.5 cortex screw,

Φ6.5 cancellous screw, Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |  |
|---------------------|------------------|------|-------|--|
| 10.14.23.05100400   | 左 L              | 5 孔  | 131mm |  |
| 10.14.23.05200400   | 右 R              | 5 孔  | 131mm |  |
| * 10.14.23.07100400 | 左 L              | 7 孔  | 167mm |  |
| 10.14.23.07200400   | 右 R              | 7 孔  | 167mm |  |
| 10.14.23.09100400   | 左 L              | 9 孔  | 203mm |  |
| 10.14.23.09200400   | 右 R              | 9 孔  | 203mm |  |
| 10.14.23.11100400   | 左 L              | 11 孔 | 239mm |  |
| 10.14.23.11200400   | 右 R              | 11 孔 | 239mm |  |

## 万向锁定股骨近端接骨板

Multi-axial proximal femur locking plate

**适用部位：**粗隆下端

Suitable for proximal femur

**配用：**

Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉

**Matching** Φ4.5 cortex screw,

Φ6.5 cancellous screw, Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.24.05100000   | 左 L              | 5 孔  | 131mm |
| 10.14.24.05200000   | 右 R              | 5 孔  | 131mm |
| * 10.14.24.07100000 | 左 L              | 7 孔  | 167mm |
| 10.14.24.07200000   | 右 R              | 7 孔  | 167mm |
| 10.14.24.09100000   | 左 L              | 9 孔  | 203mm |
| 10.14.24.09200000   | 右 R              | 9 孔  | 203mm |
| 10.14.24.11100000   | 左 L              | 11 孔 | 239mm |
| 10.14.24.11200000   | 右 R              | 11 孔 | 239mm |



## 锁定股骨远端接骨板

Distal femur locking plate

**适用部位：**股骨远端

Suitable for distal femur

**配用：**Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉

**Matching** Φ4.5 cortex screw, Φ6.5 cancellous screw, Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.26.05102400   | 左 L              | 5 孔  | 153mm |
| 10.14.26.05202400   | 右 R              | 5 孔  | 153mm |
| ★ 10.14.26.07102400 | 左 L              | 7 孔  | 189mm |
| 10.14.26.07202400   | 右 R              | 7 孔  | 189mm |
| 10.14.26.09102400   | 左 L              | 9 孔  | 225mm |
| 10.14.26.09202400   | 右 R              | 9 孔  | 225mm |
| 10.14.26.11102400   | 左 L              | 11 孔 | 261mm |
| 10.14.26.11202400   | 右 R              | 11 孔 | 261mm |

## 万向锁定股骨远端接骨板

Multi-axial distal femur locking plate

**适用部位：**股骨远端

Suitable for distal femur

**配用：**Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉

**Matching** Φ4.5 cortex screw, Φ6.5 cancellous screw, Φ5.0 locking screw

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.27.05102000   | 左 L              | 5 孔  | 153mm |
| 10.14.27.05202000   | 右 R              | 5 孔  | 153mm |
| * 10.14.27.07102000 | 左 L              | 7 孔  | 189mm |
| 10.14.27.07202000   | 右 R              | 7 孔  | 189mm |
| 10.14.27.09102000   | 左 L              | 9 孔  | 225mm |
| 10.14.27.09202000   | 右 R              | 9 孔  | 225mm |
| 10.14.27.11102000   | 左 L              | 11 孔 | 261mm |
| 10.14.27.11202000   | 右 R              | 11 孔 | 261mm |





## 锁定胫骨平台内侧接骨板 Medial tibia plateau locking plate

**适用部位：**胫骨平台内侧  
Suitable for medial tibia plateau

**配用：**Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉  
**Matching** Φ4.5 cortex screw, Φ6.5 cancellous screw, Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.28.03102000   | 左 L              | 3 孔  | 89mm  |
| 10.14.28.03202000   | 右 R              | 3 孔  | 89mm  |
| 10.14.28.05102000   | 左 L              | 5 孔  | 125mm |
| 10.14.28.05202000   | 右 R              | 5 孔  | 125mm |
| * 10.14.28.07102000 | 左 L              | 7 孔  | 161mm |
| 10.14.28.07202000   | 右 R              | 7 孔  | 161mm |
| 10.14.28.09102000   | 左 L              | 9 孔  | 197mm |
| 10.14.28.09202000   | 右 R              | 9 孔  | 197mm |
| 10.14.28.11102000   | 左 L              | 11 孔 | 233mm |
| 10.14.28.11202000   | 右 R              | 11 孔 | 233mm |

## 万向锁定胫骨平台内侧接骨板

### Multi-axial medial tibia plateau locking plate

**适用部位：**胫骨平台内侧  
Suitable for medial tibia plateau

**配用：**Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉  
**Matching** Φ4.5 cortex screw, Φ6.5 cancellous screw, Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.29.03102000   | 左 L              | 3 孔  | 89mm  |
| 10.14.29.03202000   | 右 R              | 3 孔  | 89mm  |
| 10.14.29.05102000   | 左 L              | 5 孔  | 125mm |
| 10.14.29.05202000   | 右 R              | 5 孔  | 125mm |
| * 10.14.29.07102000 | 左 L              | 7 孔  | 161mm |
| 10.14.29.07202000   | 右 R              | 7 孔  | 161mm |
| 10.14.29.09102000   | 左 L              | 9 孔  | 197mm |
| 10.14.29.09202000   | 右 R              | 9 孔  | 197mm |
| 10.14.29.11102000   | 左 L              | 11 孔 | 233mm |
| 10.14.29.11202000   | 右 R              | 11 孔 | 233mm |



## 锁定胫骨平台外侧接骨板 Lateral tibia plateau locking plate

**适用部位：**胫骨平台外侧

Suitable for lateral tibia plateau

**配用：**Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉

**Matching** Φ4.5 cortex screw, Φ6.5 cancellous screw,  
Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.30.03101000   | 左 L              | 3 孔  | 100mm |
| 10.14.30.03201000   | 右 R              | 3 孔  | 100mm |
| 10.14.30.05101000   | 左 L              | 5 孔  | 136mm |
| 10.14.30.05201000   | 右 R              | 5 孔  | 136mm |
| ★ 10.14.30.07101000 | 左 L              | 7 孔  | 172mm |
| 10.14.30.07201000   | 右 R              | 7 孔  | 172mm |
| 10.14.30.09101000   | 左 L              | 9 孔  | 208mm |
| 10.14.30.09201000   | 右 R              | 9 孔  | 208mm |
| 10.14.30.11101000   | 左 L              | 11 孔 | 244mm |
| 10.14.30.11201000   | 右 R              | 11 孔 | 244mm |

## 万向锁定胫骨平台外侧接骨板

Multi-axial lateral tibia plateau locking plate

**适用部位：**胫骨平台外侧

Suitable for lateral tibia plateau

**配用：**Φ4.5 皮质骨螺钉, Φ6.5 松质骨螺钉, Φ5.0 锁定钉

**Matching** Φ4.5 cortex screw, Φ6.5 cancellous screw, Φ5.0 locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.31.03101000   | 左 L              | 3 孔  | 100mm |
| 10.14.31.03201000   | 右 R              | 3 孔  | 100mm |
| 10.14.31.05101000   | 左 L              | 5 孔  | 136mm |
| 10.14.31.05201000   | 右 R              | 5 孔  | 136mm |
| * 10.14.31.07101000 | 左 L              | 7 孔  | 172mm |
| 10.14.31.07201000   | 右 R              | 7 孔  | 172mm |
| 10.14.31.09101000   | 左 L              | 9 孔  | 208mm |
| 10.14.31.09201000   | 右 R              | 9 孔  | 208mm |
| 10.14.31.11101000   | 左 L              | 11 孔 | 244mm |
| 10.14.31.11201000   | 右 R              | 11 孔 | 244mm |



## 锁定胫骨远端外侧接骨板

Distal lateral tibia locking plate

**适用部位：** 胫骨中下端

Suitable for tibia middle and lower end

**配用：**  $\Phi 4.5$  松质骨螺钉,  $\Phi 4.5$  皮质骨螺钉,  $\Phi 5.0$  锁定钉

**Matching**  $\Phi 4.5$  cancellous screw,  $\Phi 4.5$  cortex screw,  $\Phi 5.0$  locking screw



| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.33.05101000   | 左 L              | 5 孔  | 121mm |
| 10.14.33.05201000   | 右 R              | 5 孔  | 121mm |
| * 10.14.33.07101000 | 左 L              | 7 孔  | 157mm |
| 10.14.33.07201000   | 右 R              | 7 孔  | 157mm |
| 10.14.33.09101000   | 左 L              | 9 孔  | 193mm |
| 10.14.33.09201000   | 右 R              | 9 孔  | 193mm |
| 10.14.33.11101000   | 左 L              | 11 孔 | 229mm |
| 10.14.33.11201000   | 右 R              | 11 孔 | 229mm |

## 万向锁定胫骨远端外侧接骨板

### Multi-axial distal lateral tibia locking plate

**适用部位:** 胫骨中下端

Suitable for tibia middle and lower end

**配用:**  $\Phi 4.5$  松质骨螺钉,  $\Phi 4.5$  皮质骨螺钉,  $\Phi 5.0$  锁定钉

**Matching**  $\Phi 4.5$  cancellous screw,  $\Phi 4.5$  cortex screw,  $\Phi 5.0$  locking screw

| 订货编码 Order code.    | 规格 Specification |      |       |
|---------------------|------------------|------|-------|
| 10.14.34.05101000   | 左 L              | 5 孔  | 121mm |
| 10.14.34.05201000   | 右 R              | 5 孔  | 121mm |
| * 10.14.34.07101000 | 左 L              | 7 孔  | 157mm |
| 10.14.34.07201000   | 右 R              | 7 孔  | 157mm |
| 10.14.34.09101000   | 左 L              | 9 孔  | 193mm |
| 10.14.34.09201000   | 右 R              | 9 孔  | 193mm |
| 10.14.34.11101000   | 左 L              | 11 孔 | 229mm |
| 10.14.34.11201000   | 右 R              | 11 孔 | 229mm |





## 钛质锁定接骨螺钉

Titanium locking screw

### Φ2.0 锁定钉

Φ2.0 Locking screw

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.06.0120.007119 | 2.0 × 7               |
| 11.06.0120.009119 | 2.0 × 9               |
| 11.06.0120.011119 | 2.0 × 11              |
| 11.06.0120.013119 | 2.0 × 13              |
| 11.06.0120.015119 | 2.0 × 15              |
| 11.06.0120.017119 | 2.0 × 17              |
| 11.06.0120.019119 | 2.0 × 19              |



**配用：**

Φ1.6 医用钻头  
SW1.2 四角改锥

**Matching:**

Φ1.6 medical drill bit  
SW1.2 quadrangle  
screwdriver

### Φ2.4 锁定钉

Φ2.4 Locking screw

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.11.0524.010119 | 2.4 × 10              |
| 11.11.0524.012119 | 2.4 × 12              |
| 11.11.0524.014119 | 2.4 × 14              |
| 11.11.0524.016119 | 2.4 × 16              |
| 11.11.0524.018119 | 2.4 × 18              |
| 11.11.0524.020119 | 2.4 × 20              |
| 11.11.0524.022119 | 2.4 × 22              |
| 11.11.0524.024119 | 2.4 × 24              |
| 11.11.0524.026119 | 2.4 × 26              |



**配用：**

Φ1.9 医用钻头  
SW1.5 四角改锥

**Matching:**

Φ1.9 medical drill bit  
SW1.5 quadrangle  
screwdriver

# Φ3.0 锁定钉

## Φ3.0 Locking screw

| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.11.0530.012118 | 3.0 × 12             |
| 11.11.0530.014118 | 3.0 × 14             |
| 11.11.0530.016118 | 3.0 × 16             |
| 11.11.0530.018118 | 3.0 × 18             |
| 11.11.0530.020118 | 3.0 × 20             |
| 11.11.0530.022118 | 3.0 × 22             |
| 11.11.0530.024118 | 3.0 × 24             |
| 11.11.0530.026118 | 3.0 × 26             |
| 11.11.0530.028118 | 3.0 × 28             |
| 11.11.0530.030118 | 3.0 × 30             |
| 11.11.0530.032118 | 3.0 × 32             |
| 11.11.0530.034118 | 3.0 × 34             |
| 11.11.0530.036118 | 3.0 × 36             |
| 11.11.0530.038118 | 3.0 × 38             |
| 11.11.0530.040118 | 3.0 × 40             |
| 11.11.0530.042118 | 3.0 × 42             |
| 11.11.0530.044118 | 3.0 × 44             |
| 11.11.0530.046118 | 3.0 × 46             |



### 配用:

Φ2.3 医用钻头  
SW2.0 四角改锥

### Matching:

Φ2.3 medical drill bit  
SW2.0 quadrangle  
screwdriver



## Φ4.0 锁定钉 Φ4.0 Locking screw

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.11.0340.012112 | 4.0 × 12              |
| 11.11.0340.014112 | 4.0 × 14              |
| 11.11.0340.016112 | 4.0 × 16              |
| 11.11.0340.018112 | 4.0 × 18              |
| 11.11.0340.020112 | 4.0 × 20              |
| 11.11.0340.022112 | 4.0 × 22              |
| 11.11.0340.024112 | 4.0 × 24              |
| 11.11.0340.026112 | 4.0 × 26              |
| 11.11.0340.028112 | 4.0 × 28              |
| 11.11.0340.030112 | 4.0 × 30              |
| 11.11.0340.032112 | 4.0 × 32              |
| 11.11.0340.034112 | 4.0 × 34              |
| 11.11.0340.036112 | 4.0 × 36              |
| 11.11.0340.038112 | 4.0 × 38              |
| 11.11.0340.040112 | 4.0 × 40              |
| 11.11.0340.042112 | 4.0 × 42              |
| 11.11.0340.044112 | 4.0 × 44              |
| 11.11.0340.046112 | 4.0 × 46              |
| 11.11.0340.048112 | 4.0 × 48              |
| 11.11.0340.050112 | 4.0 × 50              |
| 11.11.0340.052112 | 4.0 × 52              |
| 11.11.0340.054112 | 4.0 × 54              |
| 11.11.0340.056112 | 4.0 × 56              |
| 11.11.0340.058112 | 4.0 × 58              |
| 11.11.0340.060112 | 4.0 × 60              |
| 11.11.0340.065112 | 4.0 × 65              |
| 11.11.0340.070112 | 4.0 × 70              |
| 11.11.0340.075112 | 4.0 × 75              |
| 11.11.0340.080112 | 4.0 × 80              |
| 11.11.0340.085112 | 4.0 × 85              |
| 11.11.0340.090112 | 4.0 × 90              |
| 11.11.0340.095112 | 4.0 × 95              |



### 配用:

Φ3.2 医用钻头  
SW2.5 六角改锥

### Matching:

Φ3.2 medical drill bit  
SW2.5 hexagon screwdriver

## Φ5.0 锁定钉 Φ5.0 Locking screw

| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.11.0350.012112 | 5.0 × 12             |
| 11.11.0350.014112 | 5.0 × 14             |
| 11.11.0350.016112 | 5.0 × 16             |
| 11.11.0350.018112 | 5.0 × 18             |
| 11.11.0350.020112 | 5.0 × 20             |
| 11.11.0350.022112 | 5.0 × 22             |
| 11.11.0350.024112 | 5.0 × 24             |
| 11.11.0350.026112 | 5.0 × 26             |
| 11.11.0350.028112 | 5.0 × 28             |
| 11.11.0350.030112 | 5.0 × 30             |
| 11.11.0350.032112 | 5.0 × 32             |
| 11.11.0350.034112 | 5.0 × 34             |
| 11.11.0350.036112 | 5.0 × 36             |
| 11.11.0350.038112 | 5.0 × 38             |
| 11.11.0350.040112 | 5.0 × 40             |
| 11.11.0350.042112 | 5.0 × 42             |
| 11.11.0350.044112 | 5.0 × 44             |
| 11.11.0350.046112 | 5.0 × 46             |
| 11.11.0350.048112 | 5.0 × 48             |
| 11.11.0350.050112 | 5.0 × 50             |
| 11.11.0350.055112 | 5.0 × 55             |
| 11.11.0350.060112 | 5.0 × 60             |
| 11.11.0350.065112 | 5.0 × 65             |
| 11.11.0350.070112 | 5.0 × 70             |
| 11.11.0350.075112 | 5.0 × 75             |
| 11.11.0350.080112 | 5.0 × 80             |
| 11.11.0350.085112 | 5.0 × 85             |
| 11.11.0350.090112 | 5.0 × 90             |
| 11.11.0350.095112 | 5.0 × 95             |
| 11.11.0350.100112 | 5.0 × 100            |



### 配用:

Φ4.2 医用钻头  
SW3.5 六角改锥

### Matching:

Φ4.2 medical drill bit  
SW3.5 hexagon  
screwdriver



## Φ6.5 空心锁定钉

### Φ6.5 Cannulated locking screw

#### 配用:

Φ2.5 螺纹导针

Φ5.0 空心医用钻头

SW4.0 空心六角改锥



#### Matching:

Φ2.5 threaded needle

Φ5.0 cannulated medical drill bit

SW4.0 cannulated hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.11.0665.070113 | 6.5 × 70              |
| 11.11.0665.075113 | 6.5 × 75              |
| 11.11.0665.080113 | 6.5 × 80              |
| 11.11.0665.085113 | 6.5 × 85              |
| 11.11.0665.090113 | 6.5 × 90              |
| 11.11.0665.095113 | 6.5 × 95              |
| 11.11.0665.100113 | 6.5 × 100             |
| 11.11.0665.105113 | 6.5 × 105             |
| 11.11.0665.110113 | 6.5 × 110             |

## 钛质皮质骨接骨螺钉

Titanium cortex screw

### Φ2.0 皮质骨螺钉

Φ2.0 Cortex screw

#### 配用:

Φ 1.6 医用钻头

SW1.2 四角改锥



#### Matching:

Φ 1.6 medical drill bit

SW1.2 quadrangle screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.07.0220.007113 | 2.0 × 7              |
| 11.07.0220.009113 | 2.0 × 9              |
| 11.07.0220.011113 | 2.0 × 11             |
| 11.07.0220.013113 | 2.0 × 13             |
| 11.07.0220.015113 | 2.0 × 15             |
| 11.07.0220.017113 | 2.0 × 17             |
| 11.07.0220.019113 | 2.0 × 19             |



## Φ2.4 皮质骨螺钉

### Φ2.4 Cortex screw

#### 配用:

Φ 1.9 医用钻头  
SW1.5 四角改锥

#### Matching:

Φ 1.9 medical drill bit  
SW1.5 quadrangle screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.12.0224.010119 | 2.4 × 10             |
| 11.12.0224.012119 | 2.4 × 12             |
| 11.12.0224.014119 | 2.4 × 14             |
| 11.12.0224.016119 | 2.4 × 16             |
| 11.12.0224.018119 | 2.4 × 18             |
| 11.12.0224.020119 | 2.4 × 20             |
| 11.12.0224.022119 | 2.4 × 22             |
| 11.12.0224.024119 | 2.4 × 24             |
| 11.12.0224.026119 | 2.4 × 26             |

## Φ3.0 皮质骨螺钉

### Φ3.0 Cortex screw

#### 配用:

Φ2.3 医用钻头  
SW2.0 四角改锥

#### Matching:

Φ2.3 medical drill bit  
SW2.0 quadrangle screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.12.0330.012118 | 3.0 × 12             |
| 11.12.0330.014118 | 3.0 × 14             |
| 11.12.0330.016118 | 3.0 × 16             |
| 11.12.0330.018118 | 3.0 × 18             |
| 11.12.0330.020118 | 3.0 × 20             |
| 11.12.0330.022118 | 3.0 × 22             |
| 11.12.0330.024118 | 3.0 × 24             |
| 11.12.0330.026118 | 3.0 × 26             |
| 11.12.0330.028118 | 3.0 × 28             |
| 11.12.0330.030118 | 3.0 × 30             |
| 11.12.0330.032118 | 3.0 × 32             |
| 11.12.0330.034118 | 3.0 × 34             |
| 11.12.0330.036118 | 3.0 × 36             |
| 11.12.0330.038118 | 3.0 × 38             |
| 11.12.0330.040118 | 3.0 × 40             |
| 11.12.0330.042118 | 3.0 × 42             |
| 11.12.0330.044118 | 3.0 × 44             |
| 11.12.0330.046118 | 3.0 × 46             |



## Φ3.5 皮质骨螺钉

### Φ3.5 Cortex screw

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.12.0335.012113 | 3.5 × 12              |
| 11.12.0335.014113 | 3.5 × 14              |
| 11.12.0335.016113 | 3.5 × 16              |
| 11.12.0335.018113 | 3.5 × 18              |
| 11.12.0335.020113 | 3.5 × 20              |
| 11.12.0335.022113 | 3.5 × 22              |
| 11.12.0335.024113 | 3.5 × 24              |
| 11.12.0335.026113 | 3.5 × 26              |
| 11.12.0335.028113 | 3.5 × 28              |
| 11.12.0335.030113 | 3.5 × 30              |
| 11.12.0335.032113 | 3.5 × 32              |
| 11.12.0335.034113 | 3.5 × 34              |
| 11.12.0335.036113 | 3.5 × 36              |
| 11.12.0335.038113 | 3.5 × 38              |
| 11.12.0335.040113 | 3.5 × 40              |
| 11.12.0335.042113 | 3.5 × 42              |
| 11.12.0335.044113 | 3.5 × 44              |
| 11.12.0335.046113 | 3.5 × 46              |
| 11.12.0335.048113 | 3.5 × 48              |
| 11.12.0335.050113 | 3.5 × 50              |
| 11.12.0335.055113 | 3.5 × 55              |
| 11.12.0335.060113 | 3.5 × 60              |
| 11.12.0335.065113 | 3.5 × 65              |



#### 配用:

Φ2.5 医用钻头  
SW2.5 六角改锥

#### Matching:

Φ2.5 medical drill bit  
SW2.5 hexagon  
screwdriver

## Φ4.5 皮质骨螺钉

### Φ4.5 Cortex screw

| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.12.0345.020113 | 4.5 × 20             |
| 11.12.0345.022113 | 4.5 × 22             |
| 11.12.0345.024113 | 4.5 × 24             |
| 11.12.0345.026113 | 4.5 × 26             |
| 11.12.0345.028113 | 4.5 × 28             |
| 11.12.0345.030113 | 4.5 × 30             |
| 11.12.0345.032113 | 4.5 × 32             |
| 11.12.0345.034113 | 4.5 × 34             |
| 11.12.0345.036113 | 4.5 × 36             |
| 11.12.0345.038113 | 4.5 × 38             |
| 11.12.0345.040113 | 4.5 × 40             |
| 11.12.0345.042113 | 4.5 × 42             |
| 11.12.0345.044113 | 4.5 × 44             |
| 11.12.0345.046113 | 4.5 × 46             |
| 11.12.0345.048113 | 4.5 × 48             |
| 11.12.0345.050113 | 4.5 × 50             |
| 11.12.0345.052113 | 4.5 × 52             |
| 11.12.0345.054113 | 4.5 × 54             |
| 11.12.0345.056113 | 4.5 × 56             |
| 11.12.0345.058113 | 4.5 × 58             |
| 11.12.0345.060113 | 4.5 × 60             |
| 11.12.0345.065113 | 4.5 × 65             |
| 11.12.0345.070113 | 4.5 × 70             |
| 11.12.0345.075113 | 4.5 × 75             |
| 11.12.0345.080113 | 4.5 × 80             |
| 11.12.0345.085113 | 4.5 × 85             |
| 11.12.0345.090113 | 4.5 × 90             |
| 11.12.0345.095113 | 4.5 × 95             |
| 11.12.0345.100113 | 4.5 × 100            |
| 11.12.0345.105113 | 4.5 × 105            |
| 11.12.0345.110113 | 4.5 × 110            |
| 11.12.0345.115113 | 4.5 × 115            |
| 11.12.0345.120113 | 4.5 × 120            |



#### 配用:

Φ3.2 医用钻头  
SW3.5 六角改锥

#### Matching:

Φ3.2 medical drill bit  
SW3.5 hexagon  
screwdriver



## 钛质松质骨接骨螺钉

Titanium cancellous screw

### Φ4.0 松质骨螺钉 (全螺纹)

Φ4.0 Cancellous screw (Fully threaded)

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.12.0440.012113 | 4.0 × 12              |
| 11.12.0440.014113 | 4.0 × 14              |
| 11.12.0440.016113 | 4.0 × 16              |
| 11.12.0440.018113 | 4.0 × 18              |
| 11.12.0440.020113 | 4.0 × 20              |
| 11.12.0440.022113 | 4.0 × 22              |
| 11.12.0440.024113 | 4.0 × 24              |
| 11.12.0440.026113 | 4.0 × 26              |
| 11.12.0440.028113 | 4.0 × 28              |
| 11.12.0440.030113 | 4.0 × 30              |
| 11.12.0440.032113 | 4.0 × 32              |
| 11.12.0440.034113 | 4.0 × 34              |
| 11.12.0440.036113 | 4.0 × 36              |
| 11.12.0440.038113 | 4.0 × 38              |
| 11.12.0440.040113 | 4.0 × 40              |
| 11.12.0440.042113 | 4.0 × 42              |
| 11.12.0440.044113 | 4.0 × 44              |
| 11.12.0440.046113 | 4.0 × 46              |
| 11.12.0440.048113 | 4.0 × 48              |
| 11.12.0440.050113 | 4.0 × 50              |
| 11.12.0440.055113 | 4.0 × 55              |
| 11.12.0440.060113 | 4.0 × 60              |
| 11.12.0440.065113 | 4.0 × 65              |
| 11.12.0440.070113 | 4.0 × 70              |
| 11.12.0440.075113 | 4.0 × 75              |
| 11.12.0440.080113 | 4.0 × 80              |
| 11.12.0440.085113 | 4.0 × 85              |
| 11.12.0440.090113 | 4.0 × 90              |
| 11.12.0440.095113 | 4.0 × 95              |



#### 配用:

Φ2.5 医用钻头  
SW2.5 六角改锥

#### Matching:

Φ2.5 medical drill bit  
SW2.5 hexagon  
screwdriver

# Φ4.0 松质骨螺钉( 半螺纹 )

## Φ4.0 Cancellous screw (Half threaded)

| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.12.0440.014213 | 4.0×14               |
| 11.12.0440.016213 | 4.0×16               |
| 11.12.0440.018213 | 4.0×18               |
| 11.12.0440.020213 | 4.0×20               |
| 11.12.0440.022213 | 4.0×22               |
| 11.12.0440.024213 | 4.0×24               |
| 11.12.0440.026213 | 4.0×26               |
| 11.12.0440.028213 | 4.0×28               |
| 11.12.0440.030213 | 4.0×30               |
| 11.12.0440.032213 | 4.0×32               |
| 11.12.0440.034213 | 4.0×34               |
| 11.12.0440.036213 | 4.0×36               |
| 11.12.0440.038213 | 4.0×38               |
| 11.12.0440.040213 | 4.0×40               |
| 11.12.0440.042213 | 4.0×42               |
| 11.12.0440.044213 | 4.0×44               |
| 11.12.0440.046213 | 4.0×46               |
| 11.12.0440.048213 | 4.0×48               |
| 11.12.0440.050213 | 4.0×50               |
| 11.12.0440.055213 | 4.0×55               |
| 11.12.0440.060213 | 4.0×60               |
| 11.12.0440.065213 | 4.0×65               |
| 11.12.0440.070213 | 4.0×70               |
| 11.12.0440.075213 | 4.0×75               |
| 11.12.0440.080213 | 4.0×80               |



### 配用:

Φ2.5 医用钻头  
SW2.5 六角改锥

### Matching:

Φ2.5 medical drill bit  
SW2.5 hexagon  
screwdriver



## Φ4.5 松质骨螺钉( 全螺纹 )

Φ4.5 Cancellous screw (Fully threaded)

### 配用:

Φ3.2 医用钻头  
SW3.5 六角改锥

### Matching:

Φ3.2 medical drill bit  
SW3.5 hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.12.0445.020113 | 4.5 × 20             |
| 11.12.0445.025113 | 4.5 × 25             |
| 11.12.0445.030113 | 4.5 × 30             |
| 11.12.0445.035113 | 4.5 × 35             |
| 11.12.0445.040113 | 4.5 × 40             |
| 11.12.0445.045113 | 4.5 × 45             |
| 11.12.0445.050113 | 4.5 × 50             |
| 11.12.0445.055113 | 4.5 × 55             |
| 11.12.0445.060113 | 4.5 × 60             |

## Φ4.5 松质骨螺钉( 半螺纹 )

### Φ4.5 Cancellous screw (Half threaded)

#### 配用:

Φ3.2 医用钻头  
SW3.5 六角改锥

#### Matching:

Φ3.2 medical drill bit  
SW3.5 hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.12.0445.025213 | 4.5 × 25             |
| 11.12.0445.030213 | 4.5 × 30             |
| 11.12.0445.035213 | 4.5 × 35             |
| 11.12.0445.040213 | 4.5 × 40             |
| 11.12.0445.045213 | 4.5 × 45             |
| 11.12.0445.050213 | 4.5 × 50             |
| 11.12.0445.055213 | 4.5 × 55             |
| 11.12.0445.060213 | 4.5 × 60             |
| 11.12.0445.065213 | 4.5 × 65             |



## Φ6.5 松质骨螺钉( 全螺纹 )

Φ6.5 Cancellous screw (Fully threaded)

### 配用:

Φ3.2 医用钻头  
SW3.5 六角改锥

### Matching:

Φ3.2 medical drill bit  
SW3.5 hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.07.0465.035113 | 6.5 × 35             |
| 11.07.0465.040113 | 6.5 × 40             |
| 11.07.0465.045113 | 6.5 × 45             |
| 11.07.0465.050113 | 6.5 × 50             |
| 11.07.0465.055113 | 6.5 × 55             |
| 11.07.0465.060113 | 6.5 × 60             |
| 11.07.0465.065113 | 6.5 × 65             |
| 11.07.0465.070113 | 6.5 × 70             |
| 11.07.0465.075113 | 6.5 × 75             |
| 11.07.0465.080113 | 6.5 × 80             |

# Φ6.5 松质骨螺钉(半螺纹)

Φ6.5 Cancellous screw (Half threaded)

## 配用:

Φ3.2 医用钻头  
SW3.5 六角改锥

## Matching:

Φ3.2 medical drill bit  
SW3.5 hexagon screwdriver

| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.07.0465.035213 | 6.5 × 35             |
| 11.07.0465.040213 | 6.5 × 40             |
| 11.07.0465.045213 | 6.5 × 45             |
| 11.07.0465.050213 | 6.5 × 50             |
| 11.07.0465.055213 | 6.5 × 55             |
| 11.07.0465.060213 | 6.5 × 60             |
| 11.07.0465.065213 | 6.5 × 65             |
| 11.07.0465.070213 | 6.5 × 70             |
| 11.07.0465.075213 | 6.5 × 75             |
| 11.07.0465.080213 | 6.5 × 80             |
| 11.07.0465.085213 | 6.5 × 85             |
| 11.07.0465.090213 | 6.5 × 90             |
| 11.07.0465.095213 | 6.5 × 95             |
| 11.07.0465.100213 | 6.5 × 100            |
| 11.07.0465.105213 | 6.5 × 105            |
| 11.07.0465.110213 | 6.5 × 110            |
| 11.07.0465.115213 | 6.5 × 115            |
| 11.07.0465.120213 | 6.5 × 120            |





## 钛质空心接骨螺钉

Titanium cannulated screw

### Φ3.5 空心接骨螺钉

Φ3.5 Cannulated screw

#### 配用:

Φ 1.2 螺纹导针

Φ 2.5 空心医用钻头

SW2.5 空心六角改锥

#### Matching:

Φ 1.2 threaded needle

Φ 2.5 cannulated medical drill bit

SW2.5 cannulated hexagon screwdriver

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0135.020214 | 3.5 × 20              |
| 11.08.0135.022214 | 3.5 × 22              |
| 11.08.0135.024214 | 3.5 × 24              |
| 11.08.0135.026214 | 3.5 × 26              |
| 11.08.0135.028214 | 3.5 × 28              |
| 11.08.0135.030214 | 3.5 × 30              |
| 11.08.0135.032214 | 3.5 × 32              |
| 11.08.0135.034214 | 3.5 × 34              |
| 11.08.0135.036214 | 3.5 × 36              |
| 11.08.0135.038214 | 3.5 × 38              |
| 11.08.0135.040214 | 3.5 × 40              |
| 11.08.0135.042214 | 3.5 × 42              |
| 11.08.0135.044214 | 3.5 × 44              |
| 11.08.0135.046214 | 3.5 × 46              |
| 11.08.0135.048214 | 3.5 × 48              |
| 11.08.0135.050214 | 3.5 × 50              |



## Φ4.0 空心接骨螺钉

### Φ4.0 Cannulated screw

#### 配用:

Φ 1.5 螺纹导针  
Φ 3.0 空心医用钻头  
SW2.5 空心四角改锥

#### Matching:

Φ 1.5 threaded needle  
Φ 3.0 cannulated medical drill bit  
SW2.5 cannulated quadrangle screwdriver



| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0240.031214 | 4.0 × 31              |
| 11.08.0240.033214 | 4.0 × 33              |
| 11.08.0240.035214 | 4.0 × 35              |
| 11.08.0240.037214 | 4.0 × 37              |
| 11.08.0240.039214 | 4.0 × 39              |
| 11.08.0240.041214 | 4.0 × 41              |
| 11.08.0240.043214 | 4.0 × 43              |
| 11.08.0240.045214 | 4.0 × 45              |
| 11.08.0240.047214 | 4.0 × 47              |



## Φ4.5 空心接骨螺钉 Φ4.5 Cannulated screw

### 配用:

Φ1.5 螺纹导针  
Φ3.0 空心医用钻头  
SW2.5 空心四角改锥

### Matching:

Φ1.5 threaded needle  
Φ3.0 cannulated medical drill bit  
SW2.5 cannulated quadrangle screwdriver



| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0245.030213 | 4.5 × 30              |
| 11.08.0245.035213 | 4.5 × 35              |
| 11.08.0245.040213 | 4.5 × 40              |
| 11.08.0245.045213 | 4.5 × 45              |
| 11.08.0245.050213 | 4.5 × 50              |
| 11.08.0245.055213 | 4.5 × 55              |
| 11.08.0245.060213 | 4.5 × 60              |

## Φ6.5 空心接骨螺钉

### Φ6.5 Cannulated screw

#### 配用:

Φ2.5 螺纹导针  
Φ5.0 空心医用钻头  
SW4.0 空心六角改锥

#### Matching:

Φ2.5 threaded needle  
Φ5.0 cannulated medical drill bit  
SW4.0 cannulated hexagon screwdriver

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0165.050213 | 6.5 × 50              |
| 11.08.0165.055213 | 6.5 × 55              |
| 11.08.0165.060213 | 6.5 × 60              |
| 11.08.0165.065213 | 6.5 × 65              |
| 11.08.0165.070213 | 6.5 × 70              |
| 11.08.0165.075213 | 6.5 × 75              |
| 11.08.0165.080213 | 6.5 × 80              |
| 11.08.0165.085213 | 6.5 × 85              |
| 11.08.0165.090213 | 6.5 × 90              |
| 11.08.0165.095213 | 6.5 × 95              |
| 11.08.0165.100213 | 6.5 × 100             |
| 11.08.0165.105213 | 6.5 × 105             |



## 钛质垫片 Titanium washer

| 订货编码 Order code.  | 规格 Specification |
|-------------------|------------------|
| 10.05.01.01000007 | Φ7               |
| 10.05.01.01003012 | Φ12              |





## Φ7.3 空心接骨螺钉

### Φ7.3 Cannulated screw

#### 配用:

Φ2.5 螺纹导针

Φ5.0 空心医用钻头

SW4.0 空心六角改锥

#### Matching:

Φ2.5 threaded needle

Φ5.0 cannulated medical drill bit

SW4.0 cannulated hexagon screwdriver

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0173.060213 | 7.3×60                |
| 11.08.0173.065213 | 7.3×65                |
| 11.08.0173.070213 | 7.3×70                |
| 11.08.0173.075213 | 7.3×75                |
| 11.08.0173.080213 | 7.3×80                |
| 11.08.0173.085213 | 7.3×85                |
| 11.08.0173.090213 | 7.3×90                |
| 11.08.0173.095213 | 7.3×95                |
| 11.08.0173.100213 | 7.3×100               |
| 11.08.0173.105213 | 7.3×105               |
| 11.08.0173.110213 | 7.3×110               |
| 11.08.0173.115213 | 7.3×115               |



## 钛质空心接骨螺钉

Titanium Cannulated Headless Compression screw

### Φ2.4 空心接骨螺钉 - 渐变加压

Φ2.4 Cannulated Headless Compression screw

#### 配用:

Φ0.9 螺纹导针

Φ1.7 空心医用钻头

SW1.5 空心六角改锥

#### Matching:

Φ0.9 threaded needle

Φ1.7 cannulated medical drill bit

SW1.5 cannulated hexagon screwdriver

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0424.008313 | 2.4 × 8               |
| 11.08.0424.009313 | 2.4 × 9               |
| 11.08.0424.010313 | 2.4 × 10              |
| 11.08.0424.011313 | 2.4 × 11              |
| 11.08.0424.012313 | 2.4 × 12              |
| 11.08.0424.013313 | 2.4 × 13              |
| 11.08.0424.014313 | 2.4 × 14              |
| 11.08.0424.016313 | 2.4 × 16              |
| 11.08.0424.018313 | 2.4 × 18              |
| 11.08.0424.020313 | 2.4 × 20              |
| 11.08.0424.022313 | 2.4 × 22              |
| 11.08.0424.024313 | 2.4 × 24              |
| 11.08.0424.026313 | 2.4 × 26              |
| 11.08.0424.028313 | 2.4 × 28              |
| 11.08.0424.030313 | 2.4 × 30              |





## Φ2.7 空心接骨螺钉 - 渐变加压

### Φ2.7 Cannulated Headless Compression screw

#### 配用:

Φ0.9 螺纹导针

Φ1.9 空心医用钻头

SW1.5 空心六角改锥

#### Matching:

Φ0.9 threaded needle

Φ1.9 cannulated medical drill bit

SW1.5 cannulated hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0427.008313 | 2.7 × 8               |
| 11.08.0427.009313 | 2.7 × 9               |
| 11.08.0427.010313 | 2.7 × 10              |
| 11.08.0427.011313 | 2.7 × 11              |
| 11.08.0427.012313 | 2.7 × 12              |
| 11.08.0427.013313 | 2.7 × 13              |
| 11.08.0427.014313 | 2.7 × 14              |
| 11.08.0427.016313 | 2.7 × 16              |
| 11.08.0427.018313 | 2.7 × 18              |
| 11.08.0427.020313 | 2.7 × 20              |
| 11.08.0427.022313 | 2.7 × 22              |
| 11.08.0427.024313 | 2.7 × 24              |
| 11.08.0427.026313 | 2.7 × 26              |
| 11.08.0427.028313 | 2.7 × 28              |
| 11.08.0427.030313 | 2.7 × 30              |

## Φ3.0 空心接骨螺钉 - 渐变加压

### Φ3.0 Cannulated Headless Compression screw

#### 配用:

Φ1.1 螺纹导针

Φ2.1 空心医用钻头

SW2.0 空心六角改锥

#### Matching:

Φ1.1 threaded needle

Φ2.1 cannulated medical drill bit

SW2.0 cannulated hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.08.0430.016313 | 3.0 × 16             |
| 11.08.0430.018313 | 3.0 × 18             |
| 11.08.0430.020313 | 3.0 × 20             |
| 11.08.0430.022313 | 3.0 × 22             |
| 11.08.0430.024313 | 3.0 × 24             |
| 11.08.0430.026313 | 3.0 × 26             |
| 11.08.0430.028313 | 3.0 × 28             |
| 11.08.0430.030313 | 3.0 × 30             |
| 11.08.0430.032313 | 3.0 × 32             |
| 11.08.0430.034313 | 3.0 × 34             |
| 11.08.0430.036313 | 3.0 × 36             |
| 11.08.0430.038313 | 3.0 × 38             |
| 11.08.0430.040313 | 3.0 × 40             |



## Φ3.5 空心接骨螺钉 - 渐变加压

### Φ3.5 Cannulated Headless Compression screw

#### 配用:

Φ1.1 螺纹导针

Φ2.5 空心医用钻头

SW2.0 空心六角改锥

#### Matching:

Φ1.1 threaded needle

Φ2.5 cannulated medical drill bit

SW2.0 cannulated hexagon screwdriver

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0435.020313 | 3.5 × 20              |
| 11.08.0435.022313 | 3.5 × 22              |
| 11.08.0435.024313 | 3.5 × 24              |
| 11.08.0435.026313 | 3.5 × 26              |
| 11.08.0435.028313 | 3.5 × 28              |
| 11.08.0435.030313 | 3.5 × 30              |
| 11.08.0435.032313 | 3.5 × 32              |
| 11.08.0435.034313 | 3.5 × 34              |
| 11.08.0435.036313 | 3.5 × 36              |
| 11.08.0435.038313 | 3.5 × 38              |
| 11.08.0435.040313 | 3.5 × 40              |
| 11.08.0435.042313 | 3.5 × 42              |
| 11.08.0435.044313 | 3.5 × 44              |
| 11.08.0435.046313 | 3.5 × 46              |
| 11.08.0435.048313 | 3.5 × 48              |
| 11.08.0435.050313 | 3.5 × 50              |



## Φ4.0 空心接骨螺钉 - 渐变加压

### Φ4.0 Cannulated Headless Compression screw

**配用:** Φ1.5 螺纹导针  
Φ2.7 空心医用钻头  
SW2.5 空心六角改锥

**Matching:** Φ1.5 threaded needle  
Φ2.7 cannulated medical drill bit  
SW2.5 cannulated hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.08.0440.020313 | 4.0 × 20             |
| 11.08.0440.022313 | 4.0 × 22             |
| 11.08.0440.024313 | 4.0 × 24             |
| 11.08.0440.026313 | 4.0 × 26             |
| 11.08.0440.028313 | 4.0 × 28             |
| 11.08.0440.030313 | 4.0 × 30             |
| 11.08.0440.032313 | 4.0 × 32             |
| 11.08.0440.034313 | 4.0 × 34             |
| 11.08.0440.036313 | 4.0 × 36             |
| 11.08.0440.038313 | 4.0 × 38             |
| 11.08.0440.040313 | 4.0 × 40             |
| 11.08.0440.042313 | 4.0 × 42             |
| 11.08.0440.044313 | 4.0 × 44             |
| 11.08.0440.046313 | 4.0 × 46             |
| 11.08.0440.048313 | 4.0 × 48             |
| 11.08.0440.050313 | 4.0 × 50             |



## 钛质空心接骨螺钉

Titanium Cannulated Compression Screw

### Φ2.4 空心接骨螺钉 - 双向加压

Φ2.4 Cannulated Compression Screw

#### 配用:

Φ0.9 螺纹导针

Φ1.7 空心医用钻头

SW1.5 空心六角改锥

#### Matching:

Φ0.9 threaded needle

Φ1.7 cannulated medical drill bit

SW1.5 cannulated hexagon screwdriver



| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0424.010314 | 2.4 × 10              |
| 11.08.0424.011314 | 2.4 × 11              |
| 11.08.0424.012314 | 2.4 × 12              |
| 11.08.0424.013314 | 2.4 × 13              |
| 11.08.0424.014314 | 2.4 × 14              |
| 11.08.0424.016314 | 2.4 × 16              |
| 11.08.0424.018314 | 2.4 × 18              |
| 11.08.0424.020314 | 2.4 × 20              |
| 11.08.0424.022314 | 2.4 × 22              |



## Φ3.0 空心接骨螺钉 - 双向加压

### Φ3.0 Cannulated Compression Screw

#### 配用:

Φ1.1 螺纹导针

Φ2.1 空心医用钻头

SW2.0 空心六角改锥

#### Matching:

Φ1.1 threaded needle

Φ2.1 cannulated medical drill bit

SW2.0 cannulated hexagon screwdriver

| 订货编码 Order code.  | 规格 (mm) Specification |
|-------------------|-----------------------|
| 11.08.0430.012314 | 3.0 × 12              |
| 11.08.0430.013314 | 3.0 × 13              |
| 11.08.0430.014314 | 3.0 × 14              |
| 11.08.0430.016314 | 3.0 × 16              |
| 11.08.0430.018314 | 3.0 × 18              |
| 11.08.0430.020314 | 3.0 × 20              |
| 11.08.0430.022314 | 3.0 × 22              |
| 11.08.0430.024314 | 3.0 × 24              |
| 11.08.0430.026314 | 3.0 × 26              |
| 11.08.0430.028314 | 3.0 × 28              |
| 11.08.0430.030314 | 3.0 × 30              |
| 11.08.0430.032314 | 3.0 × 32              |
| 11.08.0430.034314 | 3.0 × 34              |
| 11.08.0430.036314 | 3.0 × 36              |
| 11.08.0430.038314 | 3.0 × 38              |
| 11.08.0430.040314 | 3.0 × 40              |





## Φ4.0 空心接骨螺钉 - 双向加压

### Φ4.0 Cannulated Compression Screw

#### 配备:

Φ1.5 螺纹导针

Φ3.0 空心医用钻头

SW2.5 空心四角改锥

#### Matching:

Φ1.5 threaded needle

Φ3.0 cannulated medical drill bit

SW2.5 cannulated quadrangle screwdriver



| 订货编码 Order code.  | 规格 (mm)Specification |
|-------------------|----------------------|
| 11.08.0240.031311 | 4.0 × 31             |
| 11.08.0240.033311 | 4.0 × 33             |
| 11.08.0240.035311 | 4.0 × 35             |
| 11.08.0240.037311 | 4.0 × 37             |
| 11.08.0240.039311 | 4.0 × 39             |
| 11.08.0240.041311 | 4.0 × 41             |
| 11.08.0240.043311 | 4.0 × 43             |
| 11.08.0240.045311 | 4.0 × 45             |
| 11.08.0240.047311 | 4.0 × 47             |





## 江苏双羊医疗器械有限公司

Jiangsu Shuangyang Medical Instrument Co.,Ltd.

地址：江苏省张家港市杨舍镇包基双羊大厦

邮编 Postcode: 215600

ADD: Shuangyang Building, Baoji, Yangshe  
Town, Zhangjiagang City, Jiangsu Province, China

电话 Telephone: 0512-58268332

传真 FAX: 0512-58268337

E-Mail: shuangyang@jsshuangyang.com



[www.jsshuangyang.com](http://www.jsshuangyang.com)

微信号: shuangyang8